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Instructions for Use

The device you have purchased is based on medical research proving that our modern sedentary lifestyle causes damage to the human body. Insufficient walking leads to poor blood circulation and consequently, to poor health. Throughout the world doctors have been searching for a way to enhance blood circulation without resorting to potentially hazardous medications. The Lightfeet device effectively enhances the venous and arterial circulation (up to 5 times faster), thus contributing to normal body function and maintenance of good health. Lightfeet is recommended for healthy individuals who, due to work or other reasons, find themselves confined to a sedentary position for extended periods. It is also beneficial to those who suffer from poor circulation due to a physical condition (diabetes, obesity, pregnancy, recovery from surgery, etc.). Consistent use of the device will improve your blood and lymphatic circulation and prevent health problems that arise from lack of walking. You can use the Lightfeet device at home, in the office, during long train or air travel- anywhere you can sit, you can relax and enjoy "walking." Please read carefully prior to use:

- When taking the device out of the box, keep the locking knob in the closed (horizontal) position (**Figure 1**). Locking the device prevents the pedals from moving, allowing easy transport from place to place. Place the device on a level surface in front of your chair, armchair or couch. Position the device so that the Lightfeet logo is facing you.
- Unlock the pedals by rotating the locking knob to the vertical position.
- Sit down and place your feet on the pedals. Find the most comfortable sitting position and adjust the location of the device accordingly. Your feet should remain on the pedals and be able to complete the pedaling motion smoothly without effort (**Figure 2**). The pace, pedal angle, and duration of the pedaling exercise are up to you. There is no recommended limit to the "walking" activity.

Important Notices:

Safe when used as directed. Do not stand on the device! The Lightfeet device was designed for use while sitting only. Do not lock the pedals while in use. The device is not a toy. Keep it away from children. Do not use the device if damaged. Pedals should move freely without effort. Difficulty in moving pedals and/or loud noise may indicate device is damaged and you are advised to stop using the device and contact your dealer.

Cleaning Instructions:

Wipe with soapy cloth and dry in shade. Do not soak! We wish you pleasant walking and improved health!

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The LightFeet gift box contain small nylon bag holding 3 replacement components (Picture 1) which can easily be replaced in case of material fatigue, over use etc. In this manner, you can enjoy the use of LightFeet for decades of years. The spare parts bag holds short metal axis, long metal axis, and metal spring. It is quite easy to replace these spare parts, while you only need screwdrivers (flat head & Phillips head) to do it.

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Please read carefully prior to use:

Turn the LightFeet upside down and open the two screws (Picture 2) using Phillips screwdriver, once the screws have been removed, the locking (upper) plastic element will easily come out first, and then the plastic swing element (Picture 3) by sliding it softly along the blue track of the main structure, you will now have 3 plastic elements in hands (Picture 4)

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Replacing the short metal axis:

Simply slide out with a push the short axis in the middle of the swing element and replace it

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Replacing the metal spring

Open the screw at the top of the upper plastic locking element using Phillips screwdriver (Picture 5).

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Replacing the long metal axis

Remove the green button (Picture 7) that located at the side of the LightFeet using a flat screwdriver.

You only need to remove it at the side of the broken pole

Then use the flat screwdriver to push the long axis out from the inside (Picture 8) toward the round opening where the green button has been removed (Picture 9)

Push the new axis firmly inside via the Picture 9 same axis' path (Picture 10)


Abb. / Picture 9


Abb. / Picture 10