EN

PR-960 Infrared Temperature Operation Manual

Pro
Medix

Before using this infrared body temperature gun, please read this manual carefully and use it correctly. After reading, please keep this manual properly for reference at any time when necessary. (The pictures in this manual are only for reference.)

EC REP Prolinx GmbH Brehmstr. 56, 40239, Duesseldorf
Tel: 0049 2131 4051968-0 med@eulinx.eu

Important safety instructions

- 1) Please do not use the infrared body temperature gun for the purpose beyond its original design. This product is applicable to professional or domestic purposes.
- 2) Since infrared body temperature gun is not waterproof, please do not immerse it into water or other liquids. The cleaning and disinfection shall refer to the instructions in the section of Maintenance and Storage.
- 3) The infrared body temperature gun shall be placed in the clean and dry place away from the sun, with the temperature at 5°C~40°C the humidity below 85%.
- 4) Please do not use finger to touch the lens in its head.
- 5) Please do not use laser positioning beam to shoot eyes at right angles, because it can make eyes feel uncomfortable, even affect eyesight.
- 6) The obstructions, such as sweat stain, hair, hat or scarf, may cause low temperature of the measured temperature. Therefore, please use it correctly so as to make sure the correct measure results.
- 7) Please do not drop, dismantle, repair or reform the product.
- 8) Please do not keep close to strong electrostatic field or high-intensity magnetic field so as to avoid effects on the accuracy of the measured data.
- 9) In case of any problem, please do not use it continuously and contact with retailer. Please do not fix it on your own.
- 10) If the battery and product cannot be used, please handle it according to local laws and regulations.
- 11) During long period of non-use, please take out the battery so as to avoid liquid leakage due to expired battery.
- 12) Do not mix new batteries in use with semi-used batteries so as to avoid hazard or damages to product caused by the inconsistent discharge characteristics of batteries.

Warnings

- A. Please keep the infrared body temperature gun out of reach of children;
- B. Please do not put the battery into a fire;
- C. The usage of infrared body temperature gun cannot substitute for medical treatment.

Logos and signs

	BF type applied part
	Refer to operation manual
	Cautions, please refer to attached documents
	Level of protecting against liquid inlet
	Manufacturer
	Date of manufacture
	European union representative
	Serial number

	Keep dry
	Keep away from sunlight
	Product certification

1 Product introduction:

This product is handheld non contact infrared body temperature gun. It has two types of measurement modes, including body temperature and surface temperature. It can record 34 sets of body temperature data. It has functions of high and low temperature prompt, backlight display of different body temperature ranges and auto shutdown.

1.1 Intended use:

The intended use of this infrared body temperature gun is:

- Measurement of body temperature : body
- Measurement of object temperature: for instance, temperature of rice paste;
- Measurement of liquid temperature: for instance, temperature of bath water and temperature of milk.
- This product is expected to be used by medical staff in clinic or by customers in family environment.

1.2 Contraindications: None

1.3 Product features:

- 1) Accurate and fast: The measuring error can reach $\pm 0.2^{\circ}\text{C}$ and the measuring speed can reach 1 second;
- 2) Wide temperature measurement: Measuring range can reach 0.0-100°C
- 3) Multiple use: Measurement of body temperature: body temperature; Measurement of object temperature: for instance, temperature of rice paste; Measurement of liquid temperature: for instance, temperature of bath water and milk.
4. Safety and health: Non contact measurement of forehead temperature. It has a measuring distance of 5-15cm without contacting human skin (avoid cross infection)
5. Simple and useful: One key measurement, one key switching mode.
6. Acousto-optic alarm: Green, yellow and red backlights are used to indicate whether body temperature is normal, low fever or high fever, together with sound alert of "Beep".
7. Multiple sets of memory: It can store 34 sets of measured data for easy analysis and comparison.
8. Laser positioning beam: It is convenient for user to make sure measure position.
9. Correction setup: Its setup parameters can be revised so as to adapt to ethnic sets with different skin colors and human bodies with different characteristics.

1.4 Product structure

Diagram 1 Structure diagram of infrared temperature body

1. Laser: Red laser outlet
2. Sensor: This terminal aligns the measured object in case of measurement.
3. Trigger key: It is used for power-on and measurement.
4. Up key: It is used to memory query data.
5. SET key: Short press of it can be used as On/Off key of laser positioning and long press of it can be used as On/Off key of sound.
6. Down key: It is used to memory query data.
7. MODE key: Mode switching key is used to switch the mode of body temperature and surface temperature.
8. Display screen: It displays measured data.
9. Battery cover: It can be opened to replace battery.
10. Loudspeaker

1.5 Instruction of display screen


Diagram 2 Interface of display screen

1. Body temperature
2. Surface temperature
3. Measured value of temperature
4. On/Off tip of sound
5. On/Off tip of laser positioning
6. Tip of low battery power
7. Tip of measured results
8. Unit of degree Celsius
9. Unit of Fahrenheit
10. Number of memory sets
11. Memory symbol

2. Technical specifications

Working environment: Ambient temperature: 16~35°C
Relative humidity: <85%
Battery voltage: DC3V (2xAAA)
Product size: About 150mmx75mmx50mm (length x width x height)
Product weight: 120g (without batteries)
Measuring ranges: Mode of body temperature: 34~42.9°C
Mode of surface temperature: 0~100.0°C
Measuring accuracy: 0.0~ 33.9°C $\pm 2^{\circ}\text{C}$
34.0~34.9°C $\pm 0.3^{\circ}\text{C}$
35.0~42.0°C $\pm 0.2^{\circ}\text{C}$
42.1~42.9°C $\pm 0.3^{\circ}\text{C}$
43.0~100°C $\pm 2^{\circ}\text{C}$
Measuring distance: 5~15cm
Measuring time: <2 seconds
Three color backlight:
Mode of surface temperature: Green
Mode of body temperature: 34.0~37.3°C: Green
37.4~38.0°C (set value 0.1°C): Yellow
38.1 (set value) ~42.9°C: Red
Auto power off: ≤ 10 seconds
Storage and transportation: 20°C~45°C, relative humidity: <85%

3. Instrument operation and usage

3.1 Installing batteries

- 1) Push the battery cover along the direction of arrow to open battery cover;
- 2) Put 2 AAA batteries in it and pay attention to positive and negative electrode;
- 3) Press the battery cover and push it along the opposite direction of arrow to close battery cover.

Open battery cover

Close battery cover

4

3.2 Preparing for measurement

Correct usage method is the key to measuring accuracy, or measuring error may be caused. To avoid incorrect measurement, please operate according to the following tips.

- 1) In case of measuring body temperature, infrared body temperature gun is aimed at the middle of forehead (above eyebrow center) and kept vertical, with the distance of 5 to 15cm, as shown in the following picture. When pressing the "trigger key", the measured value of temperature is displayed on the display screen at once.


- 2) Before measuring, please keep the forehead free from the coverage of hair, sweat, cosmetic or hat.
- 3) When the measured person has relatively large difference in the temperature of measuring environment, it is required to stay in testing environment for more than 5 minutes. After the temperature is consistent with ambient temperature, measuring can be conducted again.
- 4) For fever patients, cold compress of forehead and other cooling measures can be applied to make low measured data so as to avoid measuring under this situation.
- 5) The ambient temperature around the measured person shall be stable. The test cannot be conducted in the places with larger airflow, including fan and air outlets of air conditioner.
- 6) When infrared body temperature gun is taken from place with larger difference in ambient temperature to be tested, test shall be conducted after infrared body temperature gun is placed in operating environment for thirty minutes.
- 7) The infrared body temperature gun cannot be used in strong sunshine.
- 8) It is recommended to measure for about three times so as to display the highest group of data.
- 9) When measuring the temperature of human forehead, please select the mode of "body temperature"; when measuring the temperature of other object, liquid and food, please select the mode of "surface temperature".
- 10) In case of low temperature of forehead caused by some reason, you can try to aim it at ear for measurement, as shown in the following picture.

5

3.3 Starting measurement

- 1) Aim at the measured target, press the "trigger key". LCD will display all numbers and characters. Since infrared body temperature gun conducts self inspection at the time of power on, this picture displays for about 2 seconds.
- 2) After the completion of display, the "beep" will be heard (under the status that the sound of instrument is opened), which indicates that power on and measurement have finished. Meanwhile, LCD displays the temperature value of the measured target and the backlight displays one color of green, yellow and red according to the corresponding settings

6

Caution


When temperature exceeds the set value of low fever, the LCD backlight is yellow; when temperature exceeds the set value of high fever, the LCD backlight is red and LCD displays a crying face, accompanying with 3 hurried "Beep" sounds in order to prompt the patient on possible fever. To make sure the accuracy of measurement, it is required to wait for at least 30 seconds after continuously measuring for 5 times.

3.4 Switching of measurement patterns

Under the status of power-on, short press of MODE key can realize cyclic switching between the mode of "body temperature" and the mode of "surface temperature". The mode of "body temperature" can be used to measure the temperature in human body, while the mode of "surface temperature" can be used to measure the temperature of the surface of object. (Factory default is the mode of "surface temperature".)

3.5 Viewing of memory values


Under the status of power-on, short press of up key or down key can realize memory query. If you want to look over earlier measured values, shortly press down key. If you want to look over recent measured values, shortly press up key. Note: 1) The bigger the numerical value of the serial number of memory is, the earlier the measured value is. The smaller the numerical value of the serial number of memory is, the more recent the measured value is. 2) If there is memory value, it displays the memorized value directly, or it displays "—", as shown in the following diagram.

- 3) The infrared body temperature gun can store 34 sets of recently measured values at most. The values beyond 34 sets are covered automatically in chronological order.

3.6 Elimination of memory values

Under the status of power-off, press and hold MODE key, while press trigger key. It displays CLR on the screen and shuts down after 2 seconds. All the memory values are eliminated.

7

3.7 Unit setup and correction setup

The default setting parameters of this product can be revised so as to adapt to populations with different skin colors or human bodies with different characteristics as well as populations under different ambient temperatures. Default settings are conducted before delivery. If not necessary, it is not recommended to revise factory default. If revision is required, please revise according to the following methods.

- 1) Unit setup F1: Under the status of power-on, press MODE key for a long time, F1 appears on the screen. Short press of the down key can realize the switching between temperature units $^{\circ}\text{C}/^{\circ}\text{F}$.
- 2) Correction setup: Under the status of power-on, press MODE key for a long time, F1 appears on the screen, followed by another press of MODE key. The screen is switched from F1 to F2. Then, press the up key to reduce 0.1°C and press the down key to add 0.1°C.

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