axion GmbH
Mollenbachstrasse 13
71229 Leonberg
Germany
www.axion.shop

Gebrauchsanweisung/Instructions for use/Mode d'emploi/
Instrucciones de uso/Istruzioni per l'uso/Instrukcja obsługi/
Gebruiksaanwijzing

AXION STIM-PRO I-3000

EMS Beckenbodentrainer/EMS Pelvic floor trainer/
Entraîneur EMS du plancher pelvien/ Entrenador EMS de suelo pélvico/
Allenatore EMS del pavimento pelvico/
Trener EMS mięśni dna miednicy/ EMS Bekkenbodentrainer

CE 2460

V 1.0_12-2025

INHALTSVERZEICHNIS / CONTENT

4 – 26	Gebrauchsanweisung
27 – 49	Instructions for use
50 – 72	Mode d'emploi
73 – 95	Instrucciones de uso
96 – 118	Istruzioni per l'uso
119 – 140	Instrukcja obsługi
141 – 163	Gebruiksaanwijzing

DE

EN

FR

ES

IT

PL

NL

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EN Instructions for use

Quick Start Guide	28
General Information	29
Important Notices	31
Intended Use	31
Contraindications	31
Side Effects	32
Warnings & Safety Instructions	33
Program Overview	35
Pelvic Floor Treatment	36
Charging the Device	40
Operating the Device	41
Additional Functions	43
Maintenance, Transport, Storage	45
Standard Delivery Scope	45
Label	46
Function and Safety Test	47
Repairing Minor Faults	47
Technical Description	48
Warranty	49
Conformity to Safety Standards	49

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



QUICK START GUIDE

GENERAL INFORMATION

The EMS device STIM-PRO I-3000 is a muscle stimulator specifically designed for safe and effective EMS muscle stimulation for incontinence treatment and pelvic floor training. It uses gentle electrical signals delivered through a probe-electrode or electrodes to stimulate and contract the pelvic floor muscles. This promotes muscle strength as well as pelvic floor muscle and bladder control, helping to reduce or eliminate urine leakage.

Functional Principle of EMS

EMS (Electrical Muscle Stimulation) is a method that uses the positive effects of low- or medium-frequency electrical currents, which are safe for humans when applied correctly. In this type of application, electrical impulses stimulate the muscles to contract more strongly, helping to strengthen them over time.

How EMS works in Pelvic Floor Training

1. Targeted Muscle Activation:

Electrical impulses are sent through electrodes to the pelvic floor muscles, causing them to contract. This allows intense training and strengthening of the pelvic floor, even when consciously controlling the muscles is difficult.

2. Deep-acting Stimulation:

EMS also reaches deeper muscle layers which are often not fully activated during conventional exercises.

What are Pelvic Floor Muscles?

The pelvic floor consists of multiple muscle layers that support key organs such as the bowel, bladder, urethra, and uterus.

Functions of the Pelvic Floor Muscles:

- ▶ Provide support for the abdominal and pelvic organs.
- ▶ Work alongside the abdominal and back muscles to stabilize and support the spine.

Types of Incontinence

Bladder weakness and incontinence are common concerns for both men and women, impacting their long-term health. Strengthening the pelvic floor muscles is a well-established method for preventing and managing incontinence and pelvic floor weakness.

1. Stress Incontinence

Stress incontinence occurs when there is involuntary urine leakage due to an increase in abdominal pressure, such as coughing, sneezing, laughing, or physical exertion. The cause is often a weakness in the pelvic floor muscles, which are unable to adequately support the urethral closure mechanism.

A common cause of stress incontinence is a defect in the closure apparatus, which usually results from untrained or overstretched pelvic floor muscles, for example after heavy childbirth. Age-related hormone deficiencies in women, particularly a decline in estrogen levels, can weaken pelvic floor muscles and contribute to incontinence. In men, hormonal changes or a decrease in testosterone levels after surgery, such as prostate removal, may also lead to bladder control issues.

2. Urge Incontinence

Urge incontinence is characterized by sudden, intense urges to urinate, so strong that the affected person cannot reach the bathroom in time. Common causes include an overactive bladder or neurological conditions that lead to uncontrolled contractions of the bladder muscle. Those affected often feel an urgent need to urinate, regardless of physical activity.

3. Mixed Incontinence

Mixed incontinence involves symptoms of both stress and urge incontinence. Individuals experience involuntary urine loss during physical exertion as well as sudden, intense urges to urinate. This combination can significantly impact quality of life and requires a personalized therapy plan that addresses both components.

IMPORTANT NOTICES

EN

- ▶ Read the instructions for use before use. Make sure that you observe all warnings and safety instructions. Ignoring the instructions may result in user injury or device damage.
- ▶ Report all serious incidents that have occurred in connection with the product to the manufacturer and the competent authority of the Member State (where the user is located).
- ▶ STIM-PRO I-3000 delivers energy to the human body by means of electrodes or probe-electrodes for the treatment of incontinence. The CE marking with identification number 2460 applies exclusively to the EMS device and the probe- electrodes. The electrodes are CE-marked as Class I medical devices.
- ▶ Please also refer to the instructions for use provided with the electrodes and probe-electrodes you are using, paying special attention to the included safety guidelines.
- ▶ **The device has no curative effect.**

INTENDED USE

STIM-PRO I-3000 is a medical device to treat symptoms of stress, urge and mixed urinary incontinence as well as faecal incontinence. It is intended for use by adults aged 18 and older of any gender who can control the device and understand the instructions for use. It is designed to be used in homecare environment.

CONTRAINDICATIONS

- ▶ Patients with an implanted electronic device, such as a pacemaker or defibrillator, or any other metallic or electronic implant may not undergo treatment without first consulting a physician.
- ▶ Patients with heart disease, epilepsy, cancer, fever or other illnesses may not undergo treatment without prior medical consultation.
- ▶ Use the electrodes only on healthy, uninjured skin or mucosa.
- ▶ If heart disease is suspected, the treatment should only be performed under medical supervision.
- ▶ Do not stimulate via the carotid sinus nerve, especially in patients with a known sensitivity to the carotid sinus reflex.
- ▶ The stimulation must not be applied on swollen, infected or inflamed areas or skin eruptions, such as phlebitis, thrombophlebitis, varicose veins, etc.

- ▶ The stimulation must not be applied on or in the vicinity of cancerous cells.
- ▶ The stimulation must not be used in cases of numbness/ sensory disturbances in the area of application.
- ▶ Do not use the device during pregnancy.
- ▶ Do not use in cases of severe local inflammation.
- ▶ Do not use if you are using a contraceptive coil.
- ▶ Do not use in cases of higher-grade hemorrhoidal disease (when using an anal probe-electrode).
- ▶ The device should not be used in cases of bladder overflow incontinence.
- ▶ In cases of fistulae or prolapsed uterus/vagina, the device must not be used.
- ▶ The device must not be used if there are infections in the vaginal or anal area.
- ▶ The electrodes and probe-electrodes should only be used by one person to avoid the risk of cross-infections.
- ▶ The device must not be used in the case of complete denervation of the pelvic floor.
- ▶ Do not use the device within three months after pelvic surgeries, Caesarean sections or vaginal birth.

SIDE EFFECTS

The side effects of the EMS device are generally mild, even with prolonged use. If disturbing effects occur, interrupt the treatment and have a physician evaluate the further use of the device.

Skin irritation and minor burns are possible side effects. If these occur, discontinue use and consult your doctor.

If muscle stimulation is too intense and prolonged, muscle pain in the sense of muscle soreness may occur. To avoid an increase in pain, shorter and less intense stimulation should be used, especially during the first treatment sessions. During intense higher-frequency stimulation (from approx. 15 Hz), persistent muscle tension may occur. If the electrodes are applied in the chest area, this can lead to disruption of respiratory activity during stimulation.

EN

WARNINGS & SAFETY INSTRUCTIONS

EN

- ▶ In case of uncertainties regarding the use of the product by laypersons, it is essential to consult with healthcare professionals for clarification. Medical questions, especially regarding contraindications, can only be answered by the treating physician or general practitioner.
- ▶ The stimulation emitted by the device may be sufficient to cause an electric shock. Electric current of this magnitude must not flow through the rib cage (thorax) or through the chest and the upper back, as it can cause cardiac arrhythmia.
- ▶ **Do not place electrodes over the heart!**
- ▶ Do not place electrodes on the front of the throat as spasm of the laryngeal and pharyngeal muscles may occur. Stimulation over the cervical cavity (throat area) may close the airway, make breathing difficult, and have negative effects on the heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or in places where the electric current can flow transcranially (through the head).
- ▶ If the therapy is ineffective or becomes uncomfortable, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ The device does not have AP/APG protection. Do not use it in an explosive atmosphere or close to flammable mixtures.
- ▶ Electrodes must not be placed over the eyes, the mouth or inside the body.
- ▶ The device is to be used only when awake.
- ▶ Do not use the device near other electrical equipment or power sources.
- ▶ If you experience pain, skin irritation, dizziness, difficulty breathing or other unusual symptoms while using the device, discontinue use and consult a physician immediately.
- ▶ The device must only be used with vaginal/anal probe-electrodes, electrodes, or accessories that are recommended by the manufacturer.
- ▶ Do not use this device while driving, operating machinery, near water, or during activities where involuntary muscle contractions may expose the user to excessive risk of injury.
- ▶ Always turn off the unit before applying to or removing accessories from the skin or mucosa.

- ▶ Keep this device out of the reach of children.
- ▶ Electronic monitoring devices (such as ECG monitors and ECG alarms) may not function properly when electrical stimulation is used.
- ▶ **The device should be cleaned with soapy water on a clean cloth after each use.**
- ▶ The device is an electronic product and must therefore not be disposed of with residual waste. Further information at: <https://axion.shop/en/pages/environment-and-disposal>
- ▶ For first-time users, muscle stimulation can be an unfamiliar sensation. It is recommended to start at a low intensity in a sitting position to familiarize yourself with the feeling before increasing the intensity.
- ▶ Please wait until your menstrual period is over, because during this time heavy use of the abdominal muscles is not recommended.
- ▶ If the battery fluid leaks and has come into contact with the skin, rinse the affected area thoroughly with plenty of clean water.
- ▶ If the battery fluid comes in contact with your eyes, avoid rubbing your eyes. Flush your eyes thoroughly with plenty of water and seek medical attention immediately.
- ▶ If you notice unusual heat, odors, or flames coming from the battery, stop using the device.
- ▶ Simultaneous connection of a patient to a high frequency surgical equipment and the stimulator may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- ▶ Do not operate the product within 1 m of shortwave or microwave therapy equipment as it may produce instability in the stimulator output.

EN

	Aim	Description	Suggested Treatment time	Frequency of use	Intensity	Time to improvement
P1 	Treat stress incontinence	The stimulation activates muscle contractions, helping to strengthen them. Frequency : 50 Hz Pulse width: 300 µs	20 minutes	Once per day	Increase until muscles contract	Improvements may be noticeable in 4 weeks.
P2 	Treat urge incontinence	The mild stimulation calms the bladder muscle, reducing involuntary contractions and preventing unexpected leakage. Frequency : 10 Hz Pulse width: 200 µs	Continuous operation (C) preset. Start with 20 min. and increase up to 45 min. if no improvement occurs.	Once per day	Choose a convenient intensity. No muscle contractions are needed.	Improvements may be noticeable after 2 weeks
P3 	Treat mixed incontinence	A combination of stress and urge programs designed to treat mixed incontinence. Frequency: 10 / 50 Hz Pulse width: 200 / 300 µs	20 minutes (non-adjustable)	Once per day	First 10 minutes: Refer to P2 (Urge) guidance Second 10 minutes: Refer to P1 (Stress) guidance	Improvements may be noticeable after 4 weeks
P4 	After-care for maintaining pelvic muscle strength	An after-care program designed to maintain and continuously train the pelvic floor muscles after initial improvement. This program is recommended for long-term muscular support and functional stability. Frequency: 35 Hz Pulse width: 250 µs	20 minutes	Twice per week	Increase until muscles contract.	Programs should be used to maintain pelvic muscle strength

PELVIC FLOOR TREATMENT

Strengthening the Pelvic Floor Muscles through Stimulation

There are various sizes and shapes of anal and vaginal probe-electrodes to accommodate anatomical differences, individual needs and comfort requirements of men and women during pelvic floor training.

Against stress incontinence, women can use either a vaginal or a anal probe-electrode. A vaginal probe-electrode is recommended for the treatment of overactive bladder and urge incontinence.

Note:

We would like to point out that it is recommended to consult a medical professional to select the appropriate anal or vaginal probe-electrode. A well-suited probe-electrode ensures effective, comfortable and safe use to strengthen and rehabilitate the pelvic floor.

The programs are for treatment with self-adhesive electrodes or vaginal/anal probe-electrodes.

- ▶ The programs involve noticeable but still comfortable contractions of the pelvic floor muscles.
- ▶ To ensure optimal stimulation, we recommend the use of contact gel with each application.
- ▶ The standard pelvic floor training session lasts 20 minutes and should be performed once per day.

Treatment with Probe-electrode

Prepare the probe-electrode

- ▶ Clean the probe-electrode before first use and after each session lukewarm soapy water and rinse it.
- ▶ It is recommended to empty your bladder before inserting the probe-electrode.
- ▶ Apply water or a water-based lubricant to the stainless steel contacts of the probe-electrode.

Connecting the probe-electrode


- ▶ Insert the plastic end of the lead wire into the output socket at the base of the device.
- ▶ Push the metal pins of the lead wires into the sockets on the probe-electrode lead.

Caution: To avoid breaking the cable, the probe-electrode must always be pulled out by the shaft!

Inserting the Probe-electrode

Warning: Ensure the device is switched off and cables are properly attached before inserting the probe-electrode.

- Find a comfortable position. The procedure should be done lying down with slightly drawn knees.

- With the probe-electrode connected to the device and adequately lubricated, gently insert it into the vagina or the anus.
- Avoid inserting the probe-electrode too deeply; do not insert beyond the start of the cable or the shaft of the probe-electrode.

Treatment with Electrode Pad


Skin preparation

The correct preparation of the skin before placement of the electrodes will prolong the life of the electrodes and reduce the risk of skin irritation.

- Determine the locations where the electrodes are to be placed.
- Clean the area with water and mild soap which does not contain any oils; do not use alcohol. Rinse the skin and dry it thoroughly.
- Trim excess body hair with scissors; do not shave.
- Apply a small amount of contact gel to the electrode placement area to ensure even stimulation.
- When removing the electrodes, always pull them off in the direction of hair growth.

Connecting and placing the electrodes

- Connect the plug connector to the electrodes before placing the electrodes on the skin.
- Plug the other end of the cable into the output socket of the device.
- Do not pull on the cables, but only on the plugs attached to the ends of the cables.
- Do not switch on the device until all electrodes or probe-electrodes and cables are properly attached.

- **For programs P1, P3, and P4:** Place the electrodes on the perineum (the area between the anus and the genitals), as shown

in Fig. 1. Adjust the intensity until you feel a slight contraction of the anus.

- ▶ **For the treatment of urge incontinence (Program P2):** Place the electrodes just above the buttocks on either side of the spine, as shown in Fig. 2.

EN


Fig. 1 Perineal placing of electrodes for both men and women for the programs P1, P3 and P4.



Fig. 2 Electrode placing above the buttocks for P2.

CHARGING THE DEVICE

Note: Please charge the battery completely before using the device.


- ▶ Connect the one end of the USB-C cable to the device and the other end to the power adapter.
- ▶ Plug the adapter into a standard power outlet.
- ▶ The device cannot be used during the charging process. The battery symbol flashes on the LCD display while the device is charging. Once charging is complete, a battery symbol with a fully filled bar will be displayed continuously.
- ▶ Remove the power plug from the socket and disconnect the USB cable from the device.

EN

OPERATING THE DEVICE

EN

- 1 Intensity level / intensity bar
- 2 Battery level
- 3 Parameter: pulse width (µs)
- 4 Parameter: frequency (Hz)
- 5 Pause symbol
- 6 Keylock symbol
- 7 Program/ Therapy Mode symbol
- 8 Treatment time / Timer bar

	Press for 1 seconds to switch on Press for 3 seconds to switch off.
	SET: Press to confirm program/ treatment time.
	Press to switch between programs and to increase treatment time and intensity
	Press to switch between programs and to decrease treatment time and intensity
	PAUSE: Press 1 second to pause KEY LOCK: Press 3 seconds for keylock

1. Switching on the Device

Warning: Do not switch on the device until all electrodes or probe-electrodes and cables are properly attached.

- ▶ Switch on the device by pressing  for 1 second.

2. Selecting preset Programs


- ▶ The number of the program is blinking.
- ▶ By using buttons  and  you can switch between programs and corresponding therapy modes (P1 P2, P3 and P4). The symbol for the specific type of incontinence is flashing.

3. Setting the Treatment Time

- ▶ Once the program is selected, press . The program stops flashing and the treatment time starts to flash in all programs except P3 (fixed time).
- ▶ Press  and  to select the treatment duration (if available). 20 Minutes are set per default except P2. The timer can be set from 5 minutes to 60 minutes or to continuous operation C in 5-minute increments. Once the desired treatment duration is selected, press  again to confirm.

4. Setting the Intensity/Current

- ▶ After confirming the treatment time intensity level is flashing.
- ▶ Intensity level is also indicated by the half-circle bar.
- ▶ The intensity can be adjusted from 0–80 mA.
- ▶ To increase or decrease the intensity, press  (to increase) or  (to decrease) repeatedly until the desired intensity level flashes on the display. Press and hold  or  to continuously increase or decrease the value.
- ▶ A flashing circle animation will appear to indicate that the program is processing.
- ▶ Press  to confirm your selection.

EN

5. Switching off the Device

- ▶ The device switches off automatically after the therapy session time has expired.
- ▶ To turn the device off manually, press and hold  for three (3) seconds. The display fades out and the device switches off.
- ▶ In an emergency, you can also disconnect the plug of the c connection cable from the device and then remove the probe-electrode/ electrodes.
- ▶ To prevent unpleasant electric shocks, never remove the probe-electrodes/electrodes while the device is still switched on.

EN

ADDITIONAL FUNCTIONS

Pause Function

At any time during stimulation, press the Pause key  to pause the program and press it again to resume stimulation. The intensity level will ramp up slowly.

Key Lock

Press and hold Pause/Key Lock Button  for 3 seconds to lock or unlock the device; the key lock icon appears on the display.

Note: The pause function and turning off the device  can also be used at any time, even when the key lock is activated.

Treatment Time

- ▶ The timer can be set between 5 minutes and 60 minutes, or to continuous operation (C), in 5-minute increments. During operation, the remaining session time is displayed and counted down in 1-second intervals until the session ends.
- ▶ The device switches off automatically once the therapy session time has expired.

Resetting the Intensity Level

- ▶ The intensity level is reset to 0 every time the device is turned off and at the end of a therapy session.
- ▶ The treatment will be stopped if the probe-electrode or electrode is disconnected or if contact with the skin or mucosa is lost. A yellow filled circle appears.
- ▶ If you change the program during the therapy session, the intensity level is reset to 0.

Battery Indicator

The battery indicator turns yellow when the charge falls to 70 %. It turns red at 30 %, signaling that the battery should be recharged soon.

EN

Treatment Time Memory

- ▶ Press  and  simultaneously for 3 seconds to enter and exit the treatment time memory when device is shut down.
- ▶ With the  or  key, you can view 1 to 99 sets of record data (program, treatment time per therapy session).
- ▶ By pressing  you can view the total number of treatments as well as the cumulative treatment time of all treatments stored in the treatment time memory.
- ▶ If no key is pressed for 10 seconds, the screen will exit this function automatically .
- ▶ Press  and  for 3 seconds to reset treatment time memory.


Fig. 1: single treatment memory


Fig. 2: overall treatment memory

MAINTENANCE, TRANSPORT, STORAGE

Note:

- ▶ The device must not be exposed to large quantities of liquids/water.
- ▶ After use, place the device back into the foam insert of the storage case for protection.
- ▶ Store the device in the storage case in a dry and cool place if you are not going to use it for a longer period of time.
- ▶ The fully charged battery has an operating time of 5 to 6 hours.
- ▶ If you find that the battery cannot be charged effectively, or the fully charged battery can only be used for a very short time, contact our service department.
- ▶ When not in use, store the device only with a fully charged battery. The device should be fully recharged after 6 months at the latest.

STANDARD DELIVERY SCOPE



Your STIM-PRO I-3000 is delivered with the following accessories. All parts can be reordered separately from our online shop at www.axion.shop if needed.

Quantity	Product name	Order no.
1 x	EMS device	1003
1 x	USB charger and USB-C charging cable	20075
1 x	Connection cable	20032
1 x	Storage case	100309
1 x	Instructions for Use	100301

Four CE certificated self-adhesive electrodes, article number 20007 (50 × 50 mm) are supplied separately with each device. The CE certificated probe-electrodes are not part of the standard scope of delivery and can be ordered separately.

LABEL

EN

STIM-PRO I-3000  REF 1003 UDI (01)4260280185705 (11) (21) www.axion.shop CE 2460    axion GmbH Mollenbachstr.13 71229 Leonberg Germany SN		 Manufacturer REF Article number CE 2460 Conformity marking with number of the notified body (device complies with the requirements of (EU)2017/745) MD Medical device UDI Unique Device Identification SN Serial number  Date of manufacture YYYY-MM-DD
 Caution! Check the instructions for use for important safety-related information such as warnings and precautions.		 Follow instructions for use  Consult instructions for use or consult electronic instructions for use
 Temperature limit for operation	Temperature limit for operation	IP22 IP protection class
 Temperature limit for storage/transportation	Temperature limit for storage/transportation	 Device type BF
 Keep away from sunlight	Keep away from sunlight	 Keep dry
 Do not dispose of with residual waste	Do not dispose of with residual waste	 2D datamatrix barcode: encoded EAN, date of manufacture and serial number

FUNCTION AND SAFETY TEST

Precaution: For safety reasons, the following precautions should be carried out once a week.

- ▶ Check the device for external damage: deformation of the housing or damage to the electrode output sockets.
- ▶ Check that no type plates or stickers are damaged.
- ▶ Check cables and electrodes for damages.
- ▶ A safety inspection and maintenance must be carried out by authorized personnel before the device is used for the first time and before operating it after a period of inactivity; at a minimum every 24 months.
- ▶ The instructions for use must always accompany the unit. If there are any problems or defects, please contact your retailer.

REPAIRING MINOR FAULTS

Note: If the device does not work properly, the instructions below may help identify possible faults. If the problem cannot be solved with any of these tips, contact the nearest service center.

- ▶ If a cable is damaged externally, replace it immediately with an undamaged one. Only use original accessories.
- ▶ No current flowing and yellow circle flashing on screen?
 - The flashing yellow circle will appear on the screen if the probe electrode or electrode is disconnected or if contact with the skin or mucosa is lost. At the same time, the device will automatically reduce the intensity to zero.
 - Check the cable connection between the device and the probe electrode/electrodes.
 - Improve the contact between the probe electrode and the mucosa, or the electrode and the skin, by applying a water-based lubricant or contact gel.

EN

TECHNICAL DESCRIPTION

EN

Channel	1 channel with adjustable intensity in 80 levels
Pulse form	Symmetrical, bi-phasic square pulse
Product specification	EMS
Maximum Output Current ($\pm 10\%$)	80 mA (peak) with a load of 500 Ohm
Maximum Output Voltage ($\pm 10\%$)	40 V (peak) with a load of 500 Ohm
Charging port	USB Type-C; 5V DC
Power source	Rechargeable lithium-ion battery 503035, 3.7 V, 500 mAh, 1.85 Wh
Battery indicator	Low battery level indication at 30 % remaining capacity.
Timer	5–60 minutes (in 5-minute increments) and continuous operation C
Display	LCD with backlighting
Size (L/W/H)	12.7 x 6.0 x 2.4 cm
Weight	92 grams (with battery)
Programs	4
Frequency / Pulse width ($\pm 5\%$)	P1: 50 Hz / 300 μ s P2: 10 Hz / 200 μ s P3: 10/50 Hz / 200/300 μ s P4: 35 Hz / 250 μ s
Applied part classification	Type BF
IP rating	IP22
Operating conditions	Temperature: 10° C to 40° C Relative humidity: 40–90% Air pressure: 700–1013 hPa Altitude in metres: 3000
Transport and storage conditions	Temperature: -10° C to 60° C Relative humidity: 30–95% Air pressure: 500–1060 hPa

WARRANTY

All devices come with a warranty of 24 months from the date of purchase. The warranty applies to the device only (except for the battery, which is covered by a limited warranty of 6 months) and includes spare parts and labor costs. The warranty does not apply to damage resulting from improper handling, the failure to follow the instructions for use, loss or dropping.

EN

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO I-3000 device is in compliance with the safety standards EN 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances - Requirements and tests.





Special precautions, information and characteristics on “Electromagnetic compatibility (EMC)”: Medical electrical devices are subject to special requirements and precautions with regard to electromagnetic compatibility.

Warnings:

- ▶ The use of accessories, transducers and wiring other than those specified or provided by axion of this appliance may result in increased ELECTROMAGNETIC RADIATION or decreased electromagnetic immunity of the appliance and result in improper operation.
- ▶ Portable RF communications equipment (radios), including accessories such as antenna cables and external antennas, should be used no closer than 30 cm (12 inches) to any part or line of the EMS DEVICE as specified by the manufacturer. Failure to do so may result in a reduction in the performance characteristics of the device.
- ▶ Use of this appliance adjacent to other appliances or stacked with other appliances should be avoided as this may result in improper operation. If use in the manner described above is nevertheless necessary, this appliance and the other appliances should be observed to ensure that they are operating correctly.

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