



SCHALLWERK

## Music+



### **Music+** (Reusable earplug with filter)

- Weight: 1.27 gram per piece
- Material: Silicon plugs; PC filter
- Sizing: Nominal Diameter 7-14 mm
- Non allergic to skin
- Declaration of Conformity can be found in:  
[www.schallwerk.eu](http://www.schallwerk.eu)

CE 2834

EN352-2: 2002 (EU) 2016/425



Silicon earplugs are soft, comfortable to fit to the ear canal. They are designed to reduce noise evenly to maintain the full spectrum of sound while protecting hearing.

#### Insertion:

1. Before inserting the earplugs, make sure they are clean.
  2. Hold the end of the earplug between your thumb and index finger.
  3. To ease insertion, use your available hand to reach over your head and gently pull the top of the ear outward and upward. This will straighten the ear canal to easily and effectively fit in the earplug.
  4. Gently push and twist the earplug until it seals comfortably in the ear canal.
- DO NOT over insert the earplug and ensure the tab of the earplug is protruding enough for you to easily remove the plug after use.
5. When inserting, use your thumb to apply enough pressure so that the earplug moves slowly into the ear.

#### Removal:

1. To remove, gently take hold of the end of the earplug with the thumb and index finger and slowly pull the earplug out of the ear.
2. It is advisable to clean the earplugs after use with a damp cloth before storing them in the protective case.

Remember! This should not be a painful procedure. Take care when inserting and removing the earplugs and contact a medical professional if you have any medical problems with your ears.

The effectiveness of the earplugs depends on a number of factors including the size and shape of the user's ears, and noise levels in which the earplugs are used.

#### Warning:

Keep away from infants and small children. This product is non-toxic but may interfere with breathing if caught in windpipe or if cord is wrapped around neck, which could lead to serious injury or death. Improper fitting or disregard of the above recommendation will reduce effectiveness and could result in hearing loss or injury. These ear-plugs should not be used where there is a risk that the connecting cord could be caught up during use; The ear-plugs must be fitted, adjusted and maintained in accordance with these instructions; The ear-plugs must be worn at all times in noisy surroundings; These ear-plugs are marked re-usable and they must be regularly inspected for serviceability. If the recommendations above are not adhered to the protection given by the ear plugs will be severely impaired. The earplugs may be affected by certain chemical substances. Further information should be sought from the manufacturer. For additional information please see [www.schallwerk.eu](http://www.schallwerk.eu)

#### Cleaning:

The earplugs are reusable and should be inspected regularly for serviceability. Keep earplugs in a cold dry place after using them. Wash earplugs with neutral soap and lukewarm water to clean and disinfect. Dry earplugs before inserting them in the ear. Use earplugs between -20°C and 35°C. Never let the ear plug beyond 55°C. Stop wearing the earplugs if there are signs of damage or wear, replace the damaged parts to extend the lifetime of the product. New eartips are available at [www.schallwerk.eu](http://www.schallwerk.eu)

MINIMUM SOUND ATTENUATION  
SNR=23dB H=22dB M=21dB L=21dB

Frequency [Hz]	Sound attenuation [dB]	Standard deviation [dB]	APV [dB]
63	24.1	4.2	19.9
125	26.6	6.8	19.8
250	26.0	5.6	20.4
500	26.7	5.9	20.8
1000	29.2	7.8	21.4
2000	29.9	7.4	22.5
4000	27.5	7.5	20.0
8000	32.5	8.0	24.5

Test Laboratory: PZT GmbH  
Address: An der Junkerei 48F, D-28389 Wilhelmshaven, Germany  
Testing date: 21-06-2017 to 25-07-2017

# CE 2834

EN 352-2: 2002 (EU) 2016/425

Certification Notified Body:  
CCQS Certification Services Limited  
Block 1 Blanchardstown Corporate Park,  
Ballycoolin Road,  
Blanchardstown, Dublin 15 D15 AKK1  
Dublin  
Ireland