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Wuxi Exanovo Medical Instrument Co., Ltd.,
No. 42 Xixin Road, Zhangjing, Xibei Town, Wuxi City,
214194 Jiangsu China

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595 AA The Hague, Netherlands

BODY PRODUCTS RELAX Pharma und Kosmetik GmbH
Alfred-Nobel-Str. 1-3, 50226 Frechen, Germany
Made in China
Rev.07|10/22

REF N13542

Instruction for Use

Sphygmomanometer-Aneroid Premium 1-tube version

Please read instruction for use before use.

The product is in compliance with EN ISO 81060-1:2012

Purpose:

The blood pressure monitor (also known as blood pressure gauge and Sphygmomanometer) can be used externally on the upper arm or wrist to measure a person's arterial pressure.

User Information:

- Sphygmomanometer are designed exclusively as non-invasive methods. Only for medical use.
- Measurement of blood pressure only on healthy skin of the upper arm or wrist.
- To be used only by specially trained staff.
- Do not inflate above 300 mmHg.
- Do not leave the cuff inflated on the arm or thigh. There is a danger of vascular congestion.
- Limit its use to less than two minutes.
- Keep out of the reach of children.

Operating Procedure:

1. Position the pressure cuff

Practice the following procedure several times to become familiar with the equipment. Before taking your blood pressure, plan to relax and rest for at least 15 minutes. This can reduce the error due to physical activity. If you are right-handed, slip the pressure cuff over your left arm and extend your left arm palm up. Whether sitting or lying, be sure that your arm rests at the same level as your heart. Turn the cuff so that the stethoscope diaphragm (sound head) is on the inside of the arm just above the elbow. This locates it over the large artery (arteria brachialis) near the hollow of the elbow. (This artery can be also found by feeling for the pulse of the artery) Grasp the end of the cuff and pull it snug, then wrap it around and over your arm, pressing together with Velcro tape to bind the cuff securely.

2. Inflate the Cuff

Feel the pulse of the artery with your fingertips to be sure the heads of the stethoscope is located directly over it. Insert the earpieces of the stethoscope. Close the air-flow valve on the bulb (turn clockwise) DO NOT OVER TIGHTEN. Inflate the cuff by repeatedly squeezing the bulb with your right gauge. When you can no longer hear the pulse beat, raise the pressure an additional 30 mmHg.

3. Slowly deflate the Cuff

Slowly open the air-flow valve by turning counter clockwise so that the pressure drops 2-4 mmHg with each beat of your heart. (This will usually mean a drop of one to two marks on the gauge per second) The rate of deflation is important for accurate reading. Remember, the pressure of the cuff has shut off all blood. Flow to the arm. So do not leave it inflated fully any longer than absolutely necessary.

4. Systolic Pressure

After opening the air-flow valve, listen carefully for a pulse beat. The moment you hear the faint rhythmic tapping or thumping sound of the pulse beat, note the reading on the gauge. This is your systolic (upper) blood pressure. It may be different at first for you to detect the sounds of the pulse beat, listen carefully and acquaint yourself with these unfamiliar sounds. You should also notice the needle on the gauge "Bounce" when your systolic pressure is reached.

5. Diastolic Pressure

Allow the pressure to continue dropping at the same rate as before 2-4 mmHg per second. Listen carefully with the stethoscope. The sounds you hear will change the phases. From The first sharp tagging or thud, they will soften to blowing or swishing sounds. Watch the falling needle. At the exact point when you can no longer hear the sounds, read the gauge. This is your diastolic (lower) blood pressure reading.

Maintenance and Repair:

Under normal condition, the device is recommended to be calibrated every two years. Calibration should also be done after repair. Calibration should be carried out at point 50 mmHg and 200 mmHg.

General Cleaning and Care:

Clean any soiled parts of the blood pressure monitor thoroughly before use every time. All the parts of the blood pressure monitor can be cleaned with a clean cotton cloth and 2-propanol (also known as isopropyl alcohol or isopropanol) alcohol solution, or water with a mild detergent. Do not let any fluid penetrate the interior of the blood pressure monitor. The sleeve cuff can be washed by hand with warm water and a mild detergent. Allow to air-dry; do not dry in any machine.

Legal Notice:

Please report serious incidents that have occurred in connection with this product to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

Technical Specifications:

Measuring range: 0-300 mmHg

Inherent error: ± 3 mmHg

Minimal grading of the manometer plate: 2 mmHg

Bubble size: 225 x 120 mm (length x width)

Air tightness: Atmospheric pressure falls as of 4 mmHg/min.

Wuxi Exanovo Medical Instrument Co., Ltd.,
No. 42 Xixin Road, Zhangjing, Xibei Town, Wuxi City,
214194 Jiangsu China

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