axion GmbH
Mollenbachstrasse 13
71229 Leonberg
Germany
www.axion.shop

CE 2460

V 1.2 - 11/2023

Gebrauchsanweisung/Instructions for use/Mode d'emploi/
Instrucciones de uso/Istruzioni per l'uso/Instrukcja obsługi

AXION STIM-PRO Comfort+

2-Kanal Nervenstimulator/2-channel nerve stimulator/
Neurostimulateur à 2 canaux/Estimulador nervioso de 2 canales/
Elettrostimolatore nervoso a 2 canali/2-kanalowy stymulator nerwów

INHALTSVERZEICHNIS / CONTENT

2 – 14	Gebrauchsanweisung
15 – 26	Instructions for use
27 – 39	Mode d'emploi
40 – 52	Instrucciones de uso
53 – 65	Istruzioni per l'uso
66 – 78	Instrukcja obsługi

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ALLGEMEINES

Was ist TENS?

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GENERAL INFORMATION

What is TENS?

TENS (transcutaneous electrical nerve stimulation) is a proven procedure in pain therapy with few side effects. It is used for symptomatic relief of acute and chronic intractable pain and helps to improve blood circulation in the treated areas.

Functional principle of TENS

Transcutaneous electrical nerve stimulation is a non-invasive, drug-free method of pain management. TENS uses tiny electrical pulses sent through the skin to the nerves to alter your perception of pain. TENS does not cure a physiological problem; it only helps control pain. TENS doesn't work for everyone, but for most patients it helps reduce or eliminate pain.

IMPORTANT NOTICES

- ▶ Read the instructions for use before use. Make sure that you observe all warnings and safety instructions. Ignoring the instructions may result in user injury or device damage.
- ▶ Report all serious incidents that have occurred in connection with the product to the manufacturer and the competent authority of the Member State (where the user is located).
- ▶ TENS unit supplies energy to the human body by means of electrode pads. Electrode pads are classified separately as Class I products. CE 2460 is only valid for the therapy device, and CE is valid for electrodes (class I products).

INTENDED USE

The purpose of the STIM-PRO Comfort+ is the symptomatic relief of acute and chronic intractable pain. The device is used and operated by the patient himself. Use is permitted only from the age of 18 years.

INDICATIONS/CONTRAINDICATIONS

Indications

TENS is used to treat the following pain patterns: Muscle, joint or bone pain, pain associated with diseases such as arthrosis or fibromyalgia, conditions such as low back pain, neck pain or tendinitis.

Contraindications

- ▶ Patients with an implanted electronic device, such as a pacemaker or defibrillator, or any other metallic or electronic implant should not undergo TENS treatment without first consulting a physician.
- ▶ Patients with heart disease, epilepsy, cancer, fever or other illnesses must not undergo TENS treatment without prior medical consultation.
- ▶ Use the electrodes only on healthy, uninjured skin.
- ▶ If heart disease is suspected, TENS treatment should be performed only under medical supervision.
- ▶ Do not stimulate via the carotid sinus nerve, especially in patients with a known sensitivity to the carotid sinus reflex.
- ▶ The stimulation must not be applied over swollen, infected or inflamed areas or skin eruptions, such as phlebitis, thrombophlebitis, varicose veins, etc.
- ▶ The stimulation must not be applied over or in the vicinity of cancer.

Side effects

The side effects of TENS units are generally mild, even with prolonged use. If disturbing effects occur, a physician must evaluate the further use of the device.

Skin irritation and minor burns are possible side effects. If these occur, discontinue use and consult your doctor.

WARNINGS & SAFETY INSTRUCTIONS

Warnings

- ▶ The stimulation emitted by the device may be sufficient to cause an electric shock. Electric current of this magnitude must not flow through the rib cage (thorax) or through the chest and the upper back, as it can cause cardiac arrhythmia. **Do not place electrodes over the heart!**
- ▶ Do not place electrodes on the front of the throat as spasm of the laryngeal and pharyngeal muscles may occur. Stimulation over the cervical cavity (throat area) may close the airway, make breathing difficult, and have negative effects on heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or in places where the electric current can flow transcerebrally (through the head).
- ▶ If TENS therapy is ineffective or becomes uncomfortable, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ The device does not have AP/APG protection. Do not use it in the presence of explosive atmosphere and flammable mixture.
- ▶ The device may only be used according to medical advice.
- ▶ Electrodes must not be placed over the eyes, the mouth or inside the body.
- ▶ The device is to be used only when awake.
- ▶ Do not use the device during pregnancy unless prescribed by your physician.
- ▶ Simultaneous connection of a patient to a high frequency surgical equipment and the stimulator may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- ▶ Do not operate the product within 1 m of shortwave or microwave therapy equipment as it may produce instability in the stimulator output.

Safety instructions

- ▶ Use only original accessories from the manufacturer with your device. (see scope of delivery)
- ▶ Do not use this device for undiagnosed pain symptoms before consulting a physician.
- ▶ Do not use this device while driving, operating machinery, near water, or during activities where involuntary muscle contractions may expose the user to excessive risk of injury.
- ▶ Always turn off the unit before applying or removing TENS accessories to the skin.
- ▶ Keep this device out of the reach of children.
- ▶ The device has no curative effect.
- ▶ Electronic monitoring devices (such as ECG monitors and ECG alarms) may not function properly when electrical stimulation is used.
- ▶ The device should be cleaned with soapy water on a clean cloth after each use.
- ▶ The device is an electronic product and must therefore not be disposed of with residual waste. Further information at: <https://axion.shop/en/environment-and-disposal/>

MAINTENANCE, TRANSPORT, STORAGE


- ▶ The device must not be exposed to large quantities of liquids/water.
- ▶ After use, place the device back into the foam insert of the storage case for protection.
- ▶ Remove the battery and store the device in the storage case in a dry and cool place if you are not going to use it for a longer period of time.

SCOPE OF DELIVERY

Each STIM-PRO Comfort+ includes the following accessories, which can be reordered at www.axion.shop.

QTY	Product name	Order no.
1x	TENS Device	1023
1x	4x Self-adhesive electrodes 50x50mm	20007
1x	2x Connection cable	20032
1x	Storage case	1023009
1x	9.0 V Block battery	00012
1x	Instructions for use	1023001

OPERATING THE DEVICE

1 On/Off Switch & Intensity Control

The device can be switched on and off by turning the intensity knobs on and off. The current intensity can be decreased or increased by turning the intensity knobs. The intensity can be adjusted continuously from 0–80 mA per channel.

2 Display

The display lights up as soon as the device is switched on.

3 Program select button


There are 8 preset TENS programs available. These can be selected by pressing the “Program” key. The program is automatically activated 3 seconds after selection.

4 Timer ⌚

The device has a timer which can be set for 15, 30, 60 minutes and continuous operation C (continue). The activity is automatically switched off when the time has elapsed.

5 Battery level indication

When the low battery indicator is shown the battery must be replaced.

TENS – PROGRAMS

Programs 01-08 have preset parameters and have proven helpful in practice for the body areas mentioned. However, all programs can also be used for body areas other than those recommended, as each patient responds individually to the stimulation. The type of treatment and the choice of program should be decided with your healthcare professional according to your pathology.

Positioning the electrodes

The electrodes can be positioned according to the painful area to be stimulated. We remind you that it is recommended to consult a health professional to test the positioning of the electrodes for the best possible pain relief.

TENS – Programs P 01 – P 08

P 01

Pulse frequency: modulated from 15 Hz
 Pulse width: 260 μ s
 Application duration: adjustable as required (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 02

Pulse frequency: 100 Hz
 Pulse width: 260 μ s
 Application duration: adjustable as required (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 03

Pulse frequency: 120 Hz
 Pulse width: 260 μ s
 Application duration: adjustable as required (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

TENS – Programs P 01 – P 08
P 04

Pulse frequency: 2 – 100 Hz
 Pulse width: 150 – 260 μ s
 Application duration: adjustable as required
 (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 05


Pulse frequency: 100 Hz
 Pulse width: 150 – 260 μ s
 Application duration: adjustable as required
 (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 06

Pulse frequency: 80< – > 7 Hz
 Pulse width: 260 μ s
 Application duration: adjustable as required
 (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 07


Pulse frequency: 120 Hz
 Pulse width: 150 – 260 μ s
 Application duration: adjustable as required
 (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

P 08

Cyclic change of programs 1 to 7.
 Pulse frequency: Cycle
 Pulse width: P 01 – P 07
 Application duration: adjustable as required
 (15, 30 or 60 min.)
 Intensity of the current: adjustable as required (0-80 mA)

LABEL

backside ▶


	Manufacturer	REF	Article Number
	Conformity marking with number of the notified body (device complies with the requirements of (EU)2017/745)		2D Datamatrix Barcode: GTIN, date of manufacture and batch number encoded
	Date of manufacture YYYY-MM-DD	SN	Serial number
	Attention! Check the instructions for use for important safety-related information such as warnings and precautions.		Follow the operating instructions
IP22	IP protection class		Temperature limit
	Type BF device		Temperature limit Storage/Transportation
	Protect from sunlight		Protect from moisture
	Do not dispose of with residual waste	MD	Medical device
	Follow the instructions for use		

BATTERY CHANGE

To ensure proper functioning, replace the battery when the battery icon is displayed. If you are not planning to use the device for a longer period of time, please remove the battery.

- ▶ Make sure that the unit is switched off.
- ▶ Open the battery compartment (slide the cover down).
- ▶ Remove the battery.
- ▶ Insert new battery (please pay attention to polarity / suitability).
- ▶ Close battery compartment.

FUNCTION AND SAFETY TEST

For safety reasons, the following precautions should be carried out once a week.

- ▶ Check the device for any external damage:
Check for deformation of the housing or damage to the electrode output sockets.
- ▶ Check that no type plates or stickers are damaged.
- ▶ Check the indication of the LEDs when switching on.
- ▶ Check cables and electrodes for damages.
- ▶ A safety check must be performed by authorized personnel at least every two years by the end of the month in which it is due.
- ▶ The Instructions for use must always be carried with the unit. If there are any problems or defects, please contact your retailer.

REPAIR MINOR FAULTS YOURSELF

If the device does not work properly, the instructions below may help identify the cause of the fault. If the problem cannot be solved with any of these tips, contact the customer service.

- ▶ Check cables and connections! Have all plugs been connected correctly?
- ▶ Does the display not switch on? Change batteries.
- ▶ If a cable is externally damaged, replace it immediately with a good one. Only use original accessories.
- ▶ No current is flowing? Check that the electrodes are correctly connected to the same channel and that the current is strong enough.

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO Comfort+ device is in compliance with the safety standards EN 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances - Requirements and tests.

WARRANTY

All devices have a warranty of 24 months from the date of purchase. The warranty applies to the stimulator only and covers both parts and labor relating thereto. The warranty does not apply to damage resulting from improper handling, the failure to follow the instructions for use, loss or dropping.

TECHNICAL DESCRIPTION

Channel	2 channels with adjustable intensity
Intensity	Adjustable 0–80 mA with a load of 500 Ohm on each channel
Pulse amplitude	Asymmetrical, bi-phasic square pulse
Output Voltage	0–40 V
Power Source	1 pc. 9,0V battery
Size	10.0 cm (L) x 6.3 cm (W) x 3.2 cm (H)
Weight	124 g incl. battery
Program P 01	Frequency: 15 Hz Pulse width: 260 µs
Program P 02	Frequency: 100 Hz Pulse width: 260 µs
Program P 03	Frequency: 120 Hz Pulse width: 260 µs
Program P 04	Frequency: 2–100 Hz Pulse width: 150–260 µs
Program P 05	Frequency: 100 Hz Pulse width: 150–260 µs
Program P 06	Frequency: 80< – > 7 Hz Pulse width: 260 µs
Program P 07	Frequency: 120 Hz Pulse width: 150–260 µs
Program P 08	Frequency: cyclic Pulse width: P 01–P 07
Low battery indicator	Low battery display. Replace the battery as soon as possible when this icon is displayed.
Operating condition	Temperature: 10°C to 40°C Relative humidity: 40 %–90 % Air pressure: 700 hPa–1060 hPa
Transport and storage condition	Temperature: -10°C to 60°C Relative humidity: 30 %–95 % Air pressure: 700 hPa–1060 hPa
Notice	All technical values include a tolerance of ± 5 %

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