TRAININGSANLEITUNG

www.u-trainer-fitness.com

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www.u-trainer-fitness.com

PREPARATION

Before training, articulate the letter „U“. The mouth corners move slightly forwards and the lips are pointed. Concentrate on this position while working with the **U-TRAINER®**!

Important: do not apply too much resistance or use too many of the enclosed O-Rings. The lips should be strong and capable of pointing when the U-Trainer® is fully contracted. Over time, the resistance and trainings level can be increased.

easy

medium

difficult

advanced

TRAINING SESSIONS

Place the winged-ends of the **U-TRAINER®** in the mouth corners and press them together until the lips are pointed (U-Position). After reaching the farthest U-position, return to the initial position. Repeat the exercises with short breaks in between.

INSTRUCTIONS FOR AN ISOMETRICAL TRAINING SESSION

Position the **U-TRAINER®**
Hold for **60 seconds**
Relax for **10 seconds**
Repeat **3 times**

INSTRUCTIONS FOR A DYNAMICAL TRAINING SESSION

Change between Start-Position and „U“-Position for **30 times**
Relax for **10 seconds**
Repeat **5 times**

Finally, after each trainings session lightly stretch the lips with the thumb and index finger from the middle to the mouth corners.

Every person has a different, individual face and therefore must individually select a trainings plan. The above exercises and variations above are suggestions.