Bedienungsanleitung / Manual

AXION STIM-X700

2-Kanal Nerven- und
Muskelstimulator

INHALTSVERZEICHNIS / CONTENT

2 – 25	Deutsche Bedienungsanleitung	DE
26 – 49	Instruction manual in English	EN
50 – 73	Manuel d'utilisation Français	FR
74 – 97	Manual de Instrucciones Español	ES
98 – 121	Istruzioni per l'uso Italiano	IT
122 – 145	Instrukcja Obsługi Polski	PL

EINLEITUNG

DE

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INTRODUCTION

EN

The primary functions of the STIM-PRO X700 are pain relief and pain treatment with the TENS mode and electrical muscle stimulation with the EMS mode. The device is used by the patient himself directly.

Explanation of TENS

There is nothing “magic” about **T**ranscutaneous **E**lectrical **N**erve **S**timulation (TENS). TENS is intended to be used to relieve pain. The TENS unit sends comfortable impulses through the skin that stimulate the nerve (or nerves) in the treatment area. In many cases, this stimulation will greatly reduce or eliminate the pain sensation the patient feels. Pain relief varies by individual patient, mode selected for therapy, and the type of pain. In many patients, the reduction or elimination of pain lasts longer than the actual period of stimulation (sometimes as much as three to four times longer). In others, pain is only modified while stimulation actually occurs. You may discuss this with your physician or therapist.

How TENS works

Transcutaneous Electrical Nerve Stimulation is a non-invasive, drugfree method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain.

Explanation of EMS

Electrical **M**uscle **S**timulation is an accepted and proven way of treating muscular injuries. EMS has also been

used successfully for many years in the field of high-performance sport for muscle regeneration and muscle development.

How EMS works

The brain sends impulses to a muscle which causes the muscle to work. EMS generates artificial stimulation of the muscle by the use of electrical pulses that stimulate the muscles located beneath the skin. The muscle is unable to distinguish between an impulse from the brain or one that comes from outside. The muscle reacts with a contraction and is trained passively. Then when the pulse ceases, the muscle relaxes and the cycle starts over again (stimulation, contraction and relaxation).

Indications

People use TENS to relieve pain for several different types of illnesses and conditions. They use it most often to treat muscle, joint, or bone problems that occur with illnesses such as osteoarthritis or fibromyalgia, or for conditions such as low back pain, neck pain, tendinitis, or bursitis. People have also used TENS to treat sudden (acute) pain, such as labor pain, and long-lasting (chronic) pain, such as cancer pain.

EMS is widely used in hospitals and sports clinics for the treatment of muscular injuries and for re-education of paralyzed muscles, to prevent atrophy in affected muscles and improve muscle tone and blood circulation. It works for:

Relaxation of muscle spasms, prevention or retardation of disuse atrophy, increasing local blood circulation, muscle re-education, immediate post-surgical stimulation of calf muscles to prevent venous thrombosis, maintaining or increasing range of motion.

IMPORTANT SAFETY INFORMATION

Read instruction manual before operation. Be sure to comply with all “CAUTIONS” and “WARNINGS” in the manual. Failure to follow instructions can cause harm to user or device.

EN

SIDE EFFECTS

TENS

Side effects of TENS machines are generally mild, even with extended use. If any effects are especially bothersome, a physician can evaluate the continued use of the device.

Skin irritation and minor burns are potential side effects. If this occurs, discontinue use and consult your physician.

EMS

Side effects from EMS are generally very mild. Patients may experience tenderness for a day or so afterward where muscles are already sore from inflammation. Shortterm redness and skin irritation may also occur. Those who use implanted medical devices and pregnant women should avoid this therapy.

Skin irritation and minor burns are potential side effects. If this occurs, discontinue use and consult your physician.

CAUTIONS & CONTRAINDICATIONS

TENS

- ▶ Only use original accessories from the manufacturer for your device.
- ▶ Do not use this device for undiagnosed pain symptoms until consulting a physician.
- ▶ Patients with an implanted electronic device, such as a cardiac pacemaker, implanted defibrillator, or any other metallic or electronic device should not undergo TENS treatment without first consulting a doctor.
- ▶ Patients with heart disease, epilepsy, cancer, fever or any other health condition should not undergo TENS treatment without first consulting a physician.
- ▶ Stimulation delivered by this device may be sufficient to cause electrocution. Electrical current of this magnitude must not flow through the thorax or across the chest because it may cause a cardiac arrhythmia. **Do not place electrodes above the heart!**
- ▶ Do not place electrodes on the front of the throat as spasm of the Laryngeal and Pharyngeal muscle may occur. Stimulation over the carotid sinus (neck region) may close the airways, make breathing difficult, and may have adverse effects on the heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or at any sites that may cause the electrical current to flow transcerebrally (through the head).
- ▶ This device should not be used while driving, operating machinery, close to water, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

- ▶ Always turn off the device before applying or removing TENS accessories to the skin.
- ▶ Isolated cases of skin irritation may occur at the site of electrode placement following long term application. If this occurs, discontinue use and consult your physician, if the irritation does not disappear.
- ▶ If TENS therapy becomes ineffective or unpleasant, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ Keep this device out of the reach of children.
- ▶ The units have no AP/APG protection. Do not use it in the presence of explosive atmosphere and flammable mixture.
- ▶ Only apply electrodes to healthy, uninjured skin.

EMS

- ▶ Only use original accessories from the manufacturer for your device.
- ▶ Safety of powered muscle stimulators for use during pregnancy has not been established.
- ▶ Caution should be used for patients with suspected or diagnosed heart problems.
- ▶ Caution should be used for patients with suspected or diagnosed epilepsy.
- ▶ Caution should be used in the presence of the following:
 - When there is a tendency to hemorrhage following acute trauma or fracture.
 - Following recent surgical procedures when muscle contraction may disrupt the healing process.
 - Over the menstruating or pregnant uterus; and
 - Over areas of the skin which lack normal sensation.

- ▶ Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium for sensitive skin, or alternate electrode placement.
- ▶ Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- ▶ Powered muscle stimulators should be kept out of the reach of children.
- ▶ Powered muscle stimulators should be used only with the leads and electrodes recommended for use by the manufacturer.
- ▶ Portable powered muscle stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

WARNINGS

- ▶ The unit should be used only under the continued supervision of a physician.
- ▶ The unit devices have no curative value.
- ▶ The unit is a symptomatic treatment and as such suppresses the sensation of pain, which would otherwise serve as a protective mechanism.
- ▶ The unit is not effective for pain of central origin. (This includes headache)
- ▶ Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when electric stimulation is in use.
- ▶ Caution should be used in applying electric stimulation to patients suspected of having heart disease. Further clinical data is needed to show there are no adverse results.
- ▶ Electrodes should not be placed over the eyes, in the mouth, or internally.
- ▶ Do not use while sleeping.
- ▶ Do not use during pregnancy unless directed by your physician.
- ▶ The long-term effects of chronic electrical stimulation are unknown.
- ▶ Stimulation should not be applied over the carotic sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- ▶ Stimulation should not be applied transcerebrally (through the head).
- ▶ Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.
- ▶ Stimulation should not be applied over, or in proximity to, cancerous lesions.

Steps to select a program (P)


Select P mode with

In program mode you can choose between 24 different programs for various areas of application. These are fully pre-programmed - just set the intensity for the connected channels and start the application!

The program starts immediately after selecting it.

The TENS programs run in continuous mode, recognizable by a "C" shown on the right side of the display. The EMS programs have a pre-defined duration of 30 minutes. The device ends the application automatically when the time expires. The programs can also be stopped before the end of the program duration and do not have to be run through completely.

To set a program, press the „MODE“ button until you reach P mode, recognizable by a P shown in the bottom-left corner of the display.

Press  or  to select the desired program (P01 – P24). You find an overview over all available programs on pages 34 to 37.

When you chose the program, please set the intensity for the connected channels.

Increase intensity 

Decrease intensity 

In P mode, the parameters are set by default and can not be changed!

Program mode (P)

All parameters are set by default. The selected program is automatically activated for all channels.

TENS – STIMULATION

The programs have proven to be helpful in most cases for the suggested applications. However, they might be used for the treatment of other injuries or discomforts as well. Every patient responds individually to stimulation current.

TENS – Programs P1 - P12

P01	Classic TENS-mode
	To be applied in most existing types of pain (acute and chronic). Gate-Control-Effect.
P02	Low frequent TENS (Burst)
	Release of endorphins. Application for example in case of radiating pain on arms/legs/feet, muscular deeper-lying pain, lumbago-ischialgia syndrome, knee arthrosis
P03	Modulated TENS (massaging effect)
	Against habituation effect when applied for longer periods. Can be applied for most types of pain (acute and chronic).
P04	Alternating frequencies 15/2Hz
	Application for example in case of stiff joints, neck pains, lumbago, menstrual problems.
P05	Alternating frequencies 80/2 Hz
	Perfect for a variety of different pain conditions.
P06	The special program for nausea, feeling of sickness, travel sickness, seasickness
	The electrodes should be placed above the acupuncture point C6.
P07	Migraine and sensible parts of the body
	Due to the reduced duration of pulse these applications are to be specifically used in cases of tension headache, face pain, neck pain, herpes zoster, migraine.

TENS – Programs P1 - P12

P08 Chronic pain

Application in case of pains in the neck, hip, shoulders, elbows, rheumatic pains, back pain (lumbago), menstrual problems and knee arthrosis.

P09 Rheumatic pain

Pains due to rheumatoid arthritis, radiating pain on arms/ legs/feet, through modulation suitable for long-term use.

P 10 Sweep-Mode 2 – 125 Hz

In cases of therapy-resistant complaints, stimulation both on sensory and motor basis, against habituation effect.

P 11 Frequency change (2 – 100 Hz) and pulse duration modulation

Applicable for most types of acute and chronic pain. With anti-habituation effect.

P 12 Frequency modulation

Application in case of pain in the back, knee, shoulders, legs, elbows and neck, rheumatic pains.

EN

EMS – STIMULATION

The following programs for muscle stimulation have a preset timer of 30 minutes. However, you can also stop or switch the programs beforehand according to your own needs. Unless otherwise stated, use an intensity so that the muscles contract fully and a relatively short duration to maximize the training effect.

EMS – Programs P13 - P24

P13 Muscle regeneration

Recovery of the muscles after training. Use moderate intensity.

P14 Endurance program

For large muscles.

P15 Warm-up program

For large muscles.

P16 Strength training

For large muscles.

P17 Sixpack training program

Symmetrical training of the abdominal muscles. Can be used for long and lateral abdominal muscles.

P18 Legs & Buttocks - Anti-Cellulite Program

Program against cellulite. Use moderate intensity.

P19 Muscle mass training

For large muscles.

P20 Massage & Relaxation Program

Pleasant stimulation current massage. Ideal for wellness and helpful against tension. Use moderate intensity.

P21 Endurance program

For small muscles.

P22 Warm-up program

For small muscles.

P23 Strength training

For small muscles.

P24 Muscle mass training

For small muscles.

EN

What are small and large muscle groups?

The muscles can be divided into small and large muscle groups. Here you will find some examples for a better understanding:

Large muscle groups:

- ▶ thigh musculature
- ▶ buttock musculature
- ▶ abdominal musculature
- ▶ back musculature
- ▶ chest musculature

Small muscle groups:


- ▶ calf musculature
- ▶ biceps and triceps
- ▶ shoulder musculature
- ▶ forearms
- ▶ fingers and toes

OPERATING THE DEVICE


Open the front lid of the device to reach the controls. The device can be switched on increasing the intensity with the regulators on top of the machine.

EN


ON / OFF, Regulating the intensity

The intensity can be adjusted stepless from 0-100mA by turning the intensity regulators on top of the device.


Connect the electrodes with the connecting cables. The device must be turned off when connecting the electrodes. The electrodes must be firmly attached to the skin.

MODE button

There are five manual TENS modes available: **(B)** Burst, **(N)** Normal, **(M)** / **(SD1)** / **(SD2)** Modulated, as well as two EMS modes; **(S)** Synchronous, **(A)** Asynchronous, and the **(P)** Program mode. The mode can be changed by pressing the "MODE" button. In the program mode P there are 24 pre-defined programs with fixed settings available.


Setting the parameters

In the manual modes, the parameters can be selected

by pressing the „SET“ button. The keys  and  can be used to change the values of the selected parameter.


Set frequency

The frequency is adjustable from 1 Hz to 150 Hz. The current value is shown on the LCD display. The frequency determines how many electrical impulses are transmitted through the skin every second. If no instructions regarding the frequency are given in therapy, set the control to 70 – 120 Hz.


Set pulse width

The pulse width can be adjusted from 50 μs to 300 μs . It regulates the time for which the electrical signal is applied. The current value is shown on the LCD display.

Timer ⌚


The device comes with a timer that can be adjusted from 1 to 60 minutes. It can also be set to continuous mode. The application ends automatically when the timer runs out.

Set ramp time

This setting is only available in EMS mode.


The ramp time controls the time of output current that increase from 0 to the setting value, and from the setting value back to 0. The ramp time is adjustable from 1 to 8 seconds.

Set contraction time (ON-Time)

This setting is only available in EMS mode.

The ON time controls the contraction time. By pressing the (SET) control, the contraction time is adjustable from 2 – 90 seconds. This time in-


cludes the ramp time twice, once for the ramp up and once for the ramp down

**Example: 10s ON time & 3s ramp time =
3s ramp up, 4s full intensity, 3s ramp down**

Set pause time (OFF-Time)


This setting is only available in EMS mode. The Off Time controls the time of relaxation between two cycles. By pressing the  control, the relaxation time can be adjusted. Both channels' stimulation is cycled on and off by the contraction and relaxation settings. The range is adjustable from 0 to 90 seconds. In Alternate mode, the OFF time should be equal or more than the ON time.

OFF TIME ≥ ON TIME

Program mode

The programs are chosen by pressing  or  in program mode. The programs run with pre-defined settings that can not be altered (except in intensity).

ABOUT THE DEVICE

Front

EN


Back

CHANGING BATTERIES

The battery must be replaced if the performance decreases in order to guarantee the function. Remove the battery if the device is not likely to be used for some time.

- ▶ Make sure the power is off.
- ▶ Open the battery compartment by sliding down and removing the front cover.
- ▶ Remove the battery.
- ▶ Insert the new battery. Please check that the polarity is correct!
- ▶ Close the battery compartment. Attach the cover to the bearings at the bottom of the device and slide upwards.

MAINTENANCE, TRANSPORTATION AND STORAGE OF THE DEVICE

EN

- ▶ Non-flammable cleaning solution (containing 70% of alcohol) is suitable for cleaning the device.
- ▶ Stains and spots can be removed with a cleaning agent.
- ▶ Do not submerge the device in liquids or expose it to large amounts of water.
- ▶ Return the device to the carrying box with sponge foam to ensure that the unit is well-protected before transportation.
- ▶ If the device is not to be used for a long period of time, remove the batteries and put it back into the carrying box and keep it in a cool, dry place.
- ▶ The packed TENS device should be stored and transported under the temperature range of -20°C ~ $+60^{\circ}\text{C}$, relative humidity 20% ~ 95%, atmosphere pressure 500 hPa ~ 1060 hPa.

SAFETY CHECKS

For safety reasons, review the following checklist once a week.

- ▶ Check the device for deformation of the housing or damage to the output sockets.
- ▶ Make sure that the descriptions and labels are not damaged.
- ▶ Check that the LED is on when the device is turned on.
- ▶ Check the cables and electrodes for damage.

- ▶ The device must be subjected to safety checks and maintenance by authorized technicians before use and each re-use, but at least every 24 months.
- ▶ Please consult your distributor if there are any problems with the device and accessories.
- ▶ The manual must always be carried with the device.

MALFUNCTIONS

Should any malfunctions occur while using your device, check the following points. If none of these can solve the problem, please contact the customer service.

- ▶ Check if the controls are correctly adjusted. Check if the values suggested by the physician have been adjusted.
- ▶ Check if the cables are correctly connected to the device.
- ▶ The display does not turn on? Replace the batteries.
- ▶ Check the cables for possible damage. The device may only be used when subjected to a load.
- ▶ There is no stimulation current? Check if the electrodes are correctly connected to the same channel and if the intensity is high enough.

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO X700 device is in compliance with the safety standards EN60601-1:2006/A1:2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN60601-1-2:2015 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances - Requirements and tests

TECHNICAL DESCRIPTION

01	Channel	2 channels – with adjustable intensity
02	Intensity	Adjustable, 0 – 100 mA (in steps of 1 mA) with 500 ohm load on each channel
03	Pulse amplitude	Asymmetrical, Bi-phasic square pulse
04	Output voltage	0 – 50 V
05	Power source	1 pc. 9V battery
06	Size	10.2cm (L) x 6.2cm (B) x 2.8cm (H)
07	Weight	150 grams (batteries included)
08	Pulse rate	Adjustable, from 1 – 150 Hz, 1 Hz/step
09	Pulse width	Adjustable, from 50 – 300 microseconds, 10 μ /step
10	ON time	Adjustable, from 2 – 90 seconds, 1 Sec./step
11	OFF time	Adjustable, from 0 – 90 seconds, 1 Sec./step
12	Ramp time	Adjustable, from 1 – 8 seconds, 1 Sec./step
13	Mode	5 TENS modes: B (Burst) N (Normal), M (Modulation), SD1 (40% modulation), SD2 (70% modulation) – 2 EMS modes: S (Synchronous), A (Alternate) – P (Program Mode) – 24 TENS programs (set by default)
14	TENS Burst mode B	Burst rate: Adjustable, 0.5 – 5 Hz – Pulse Width adjustable, 50 – 300 μ s – Fixed internal frequency = 100 Hz
15	TENS Normal mode N	The pulse rate and pulse width are adjustable. It generates continuous stimulation based on the setting value
16	TENS Modulation mode M	Modulation mode is a combination of pulse rate and pulse width modulation. The pulse rate and width are automatically varied in a cycle pattern. The pulse width is decreased by 50% from its original setting in 0.5 second, then the pulse rate is decreased by 50% from its original setting in 0.5 second. Total cycle time is 1second. In this mode, pulse rate (1 – 150 Hz) and pulse width (50 – 300 μ s) are fully adjustable.

17	TENS SD 1 mode	The SD1 (Strength-Duration) mode consists of automatic modulation intensity and pulse width in 40% range. The intensity is always increasing while the pulse width is decreasing and vice versa. The intensity is decreased by 40% while the pulse width is increased by 40% in 5 seconds. In the next 5 seconds, the intensity is increased by 40% while the pulse width is decreased by 40%. Total cycle time is 10 seconds. Pulse rate (1 – 150 Hz) and pulse width (50 – 300 μ s) are fully adjustable.
18	TENS SD 2 mode	The SD2 (Strength Duration) mode consists of automatic modulation intensity and pulse width in 70% range. The intensity is always increasing while the pulse width is decreasing and vice versa. The intensity is decreased by 70% while the pulse width is increased by 70% in 4 seconds. In the next 4 seconds, the intensity is increased by 70% while the pulse width is decreased by 70%. Total cycle time is 10 seconds. Pulse rate (1 – 150 Hz) and pulse width (50 – 300 μ s) are fully adjustable.
19	EMS S Synchronous mode	Stimulation of all channels occurs synchronously. The "ON" time including "Ramp Up" and "Ramp Down" time. Therefore, the setting of ON Time should be no less than two times of the "Ramp" time in this mode.
20	EMS A Asynchronous mode	Alternating stimulation. Used for the training of both agonist and antagonist. The duration of one cycle is the duration of the ON time, including the ramp time. Please note that the ON time can not be shorter than twice the ramp time.
21	TENS/EMS P program mode	See section Program mode (P) – pages 34 – 37

22	Operating condition	Temperature: 0°C - 40°C Relative humidity: 30 % - 75 % Air pressure: 700hPa - 1060hPa
23	Notice	All technical values include a tolerance of +/- 5%

EN

SCOPE OF DELIVERY


Each STIM-PRO X700 is fully equipped with the following standard equipment.

- ▶ **1 piece** **TENS EMS combo device**
- ▶ **4 pieces** **Self-adhesive electrodes**
- ▶ **1 pieces** **9V battery**
- ▶ **1 piece** **Storage case**
- ▶ **2 pieces** **Electrode cables**
- ▶ **1 piece** **Instructions manual**

WARRANTY

All TENS models carry a warranty of 24 months from the date of delivery. The warranty applies to the stimulator only and covers both parts and labor relating thereto. The warranty does not apply to damage resulting from improper handling, the failure to follow the operating instructions, loss or dropping.

Manufacturer:
axion GmbH
Mollenbachstr. 13
71229 Leonberg, Germany
www.axion.shop/en

1.  Advice for operation
2.  Degree of electrical protection Type BF
3.  Do not insert the plug into the power supply socket of 230V
4.  IP protection class
5.  Do not dispose of via residual waste
6.  Increment
7.  Decrement
8.  Read the instruction manual
9.  Direct current (DC)
10.  Manufacturer
11.  Serial number

Manual Version: 1.2

Release Date of the Manual: December 2018

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