Personen-Körper Analysewaage

Body analysis scale

PRODUKTNUMMER: 10106475

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INSTRUCTIONS

General

ATTENTION! Observe the safety and assembly instructions in order to avoid the risk of injury or damaging the product.

Only for domestic use and not for commercial use!

Important: Read these instructions carefully and thoroughly. Keep these instructions to read through at a later date. If you give the product to anyone else at some point in the future, please ensure you also pass on this manual.

Keep small parts out of the reach of children. Pay attention especially when unpacking it. Keep plastic bags and other packaging away from children reach. Danger of Suffocation!

Check the delivery for completeness. Later complaints cannot be accepted.

Check all elements and parts for damages. Despite careful checks since may happen that even the best product takes damage during the transport. In this case, do not build up of your product. Defective parts can endanger your health.

Never modify the product! Changes will void the warranty and the product may be unsafe or even dangerous.

Intended Use

The scales is designed for the weighing of people in the home. It is not suitable for use in clinics, medical practices, etc.. The maximum carrying capacity of 180 kg must not be exceeded.

Use the product only for its designed purpose of use. Manufacturer will not assume any responsibility in case of damages cause of an unappropriated use. Any modification applied by customer on the original product will have a negative impact on its safety, cause danger and warranty will be cancelled.

Safety instructions

This device is not intended for use by persons (adults and children) with physical, sensory or mental impairments or lacking experience and/or knowledge of use unless they are supervised by a person responsible for their safety or have received instructions from that person on how to use the applicator.

Children should be supervised to make sure they do not play with the device.

Do not submerge the scales in water or other liquids.

Place the scales on a stable, even surface.

Do not overload the scales. Pay attention to the "Technical Data".

Protect the scales from sharp jolts and vibrations.

Do not place the scales in the vicinity of devices that emit electromagnetic radiation, such as mobile phones, wireless devices, and the like. This can lead to incorrect display indicators or malfunctions.

Do not place any objects on the scales when you are – for example – storing them in a cupboard. The sensitive measuring electronics in the scales could be damaged by heavy objects.

Do not attempt to repair the device yourself in the event of a malfunction. Repairs should only be carried out by authorized service centers.

Do not use the appliance in a location where extreme temperatures or humidity may occur. Keep the scales away from liquids.

Risk of injury

1. Do not step onto the outer edge of any one side of the scales. Risk of tipping!
2. Do not step onto the scales with wet feet or when its surface is damp. Risk of slipping!
3. The weighing surface is made of glass. If you fall down, there is a risk of the glass breaking and thus a risk of injury from glass splinters.

Suitability

Do not use these scales if you:

- are pregnant,
- have a cardiac pacemaker,,
- are a dialysis patient,
- suffer from oedema,
- have implants, e.g. artificial limbs, metal plates, etc.

If in doubt, consult your GP.

BATTERY SAFETY INSTRUCTIONS

- Do not disassemble batteries!
 - Remove discharged batteries from the device immediately!
 - Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of this parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
 - If a battery has been swallowed seek medical attention immediately!
 - Insert the batteries correctly, observing the polarity!
 - Always keep the battery compartment well closed!
 - Remove the batteries from the device if it is not going to be used for an extended period!
 - Keep batteries out of children's reach!
 - Do not attempt to recharge these batteries!
- There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
 - Do not throw into a fire! There is a danger of explosion!
 - Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!

SPECIFICATIONS

Model: 10106475

- 2 Batteries 1.5 V Type "AAA" / "LR03"
- Mass alignment: from 75 to 225 cm (2'5" – 7'4." Ft:in)
- Minimum weight: from 2.5 kg (5.5 lb/0: 4 st)
- Maximum weight: 180 kg
- Weight measurement: 100g increments
- Body fat measurement unit: 0,1 %
- Body fat measuring range: 4,0 % - 50,0 %
- Body water measuring range: 35,0 % - 85,0 %
- Muscle mass measuring ranges: 10,0 % - 50,0 %

OPERATION

Use

The device has a very sensitive sensor. If you want to adjust the settings, simply touch the buttons gently with your finger.

If the device does not respond to touch, clean the glass surface with a soft, dry cloth.

NOTE:

The scales are factory set with the units "kg" and "cm".

To change the settings to "lb" or "st" as well as "Ft:in", proceed as follows:

- Switch the scales on by applying gentle pressure on the glass surface with your foot or toe.
- Press the down button , until the required unit of weight is displayed. The unit for the height also automatically changes.

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

Location of controls

- 1 Contact Surfaces
- 2 Display
- 3 DOWN-button
- 4 UP-button
- 5 SET-button

Display Symbol	Meaning
kg/ lb/ st	Body Weight
USER 0-9	User memory space
 	Male / Female
cm (Ft:in)	Body height
Age	Age
FAT	Body fat display
TBW	Body water display
MUS	Muscle mass display
BONE	Bone weight display
BMI	Body fat index
KCAL	Body energy

Information on Measuring Body Fat

WARNING:

Do not use the scales if you use a pacemaker or another medical implant with electronic components. The function of the implant may be impaired by the measurement.

NOTE:

The device can only function correctly if you stand on the contact surfaces with bare, dry feet! The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body. The signal is generally safe! However, please read our special safety instructions!

The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data (age, gender, body weight).

There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values. Fat is a factor that puts serious strain on the body, in particular the cardiovascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values. Fat is a factor that puts serious strain on the body, in particular the cardiovascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

Saving Personal Data

NOTE:

- The scale has default values that you can change, as described below.
- Carry out the following steps in sequence and within a few seconds. Otherwise, the device will automatically cancel input mode.
- Press the button SET. The memory space will flash on the display (USER). Confirm the setting by pressing SET or select another memory space by pressing the UP / DOWN button. Then press SET to go to the next setting.
- The symbol for "male" or "female" will flash. Confirm the setting by pressing SET or use the UP / DOWN button to choose the symbol for your gender. Then press SET to go to the next setting.
- The height setting will flash. Select your height by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The setting AGE will flash. Select your age by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET
- The default symbol for the "un-athletic" type flashes. Which type do you consider yourself to be? Select your type using the UP or DOWN button: "un-athletic", "average" or "athletic". Confirm your selection with SET.
- If only the stored position and the display 0.0 kg (0.0 lb) is still shown, input mode for this stored position is completed.

You now have the following options:

- You can change the units for weight and type. Press the UP button, until the desired unit of weight is displayed.
- The unit for the height automatically changes.
- Press the SET button to see the selected settings again.
- Allocate further memory spaces. There are 10 spaces available in total. Press SET and change the space by pressing the UP or DOWN button. Proceed further as described above.
- If you want to make changes, press SET (USER flashes). Select the stored position using the UP or DOWN button. Proceed further as described above.
- You can immediately start measuring.

Measurement and Analysis

Record your personal data in the memory.

NOTE:

- The device can only function correctly if you stand on the device with bare, dry feet!
- If necessary wipe the contact surfaces to improve results.
- Due to the significant fluctuation of body fat, it is difficult to attain a standard value. On average the optimal testing time is between 6 p.m. and 8 p.m., since the body fat remains very consistent during this period of time.
- Press the SET button. The preset stored position data appears on the display (USER flashes). Use the UP / DOWN button to select the stored position at which your data is saved.
- Wait for the USER display and 0.0 kg (0.0 lb) to stop flashing.
- Stand on the contact surfaces (1) of the scales with both feet.
- Firstly your weight will flash on the display.
- After a few seconds, the following results will be shown:
 1. - The body fat content in %
- The body water content in %
- The results analysis in the form of a bar
 2. - The muscle mass content in %
- The bone mass content in %
- The calorie requirement per day
 3. - The preset gender
- The preset type
- The preset height
- The preset age (AGE)
- The weight gain (INCREASE) / the weight loss (DECREASE) in kg (lb) in comparison to the previous measurement
- Afterwards the determined data will be shown repeatedly. Then the scale will automatically turn off.

Guideline values for body fat and body water content as a percentage (approximate figures)

Age (Years)	Women		Men		
	Body fat	Body water	Body fat	Body water	
<= 30	4,0 - 20,5	66,0 - 54,7	4,0 - 15,5	66,0 - 58,1	Underweight
	20,6 - 25,0	54,6 - 51,6	15,6 - 20,0	58,0 - 55,0	Normal
	25,1 - 30,5	51,5 - 47,8	20,1 - 24,5	54,9 - 51,9	Overweight
	30,6 - 50,0	47,7 - 34,4	24,6 - 50,0	51,8 - 34,4	Obese
> 30	4,0 - 25,0	66,0 - 51,6	4,0 - 19,5	66,0 - 55,4	Underweight
	25,1 - 30,0	51,5 - 48,1	19,6 - 24,0	55,3 - 52,3	Normal
	30,1 - 35,0	48,0 - 44,7	24,1 - 28,5	52,2 - 49,2	Overweight
	35,1 - 50,0	44,6 - 34,4	28,6 - 50,0	49,1 - 34,4	Obese

Display error Messages

If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.

Display	Meaning
Lo	Battery is low
Err2	Personal wear the shoes to measuring the body fat
EEE	The weight lies above the allowed maximum weight
ErrH	Fat content over high
ErrL	Fat content lower

Cleaning

- Do not immerse the balance in water: the electronic components can be damaged.
- Do not use a wire brush or any abrasive items.
- Do not use any acidic or abrasive detergents.
- Clean the device with a slightly moist cloth. Make sure that no humidity or water seeps into the device!
- Use a commercially sold detergent or soap for cleaning.
- Dry the device well after cleaning with a dry and soft cloth.

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Disposal:

Once item end-of-life, please dispose valuable raw materials for a proper recycling. If you are not sure how to handle it in a correct way please contact your local disposal or recycling centers for advice.

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This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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