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## LIST OF CONTENTS

Items Supplied	9
At A Glance	9
Symbols	10
Signal Words	10
Intended Use	11
When May the Device Not Be Used?	11
Safety Notices	11
Charging Up the Battery	14
Use	14
Cleaning and Storage	15
Troubleshooting	15
Disposal	15
Technical Data	16

## ITEMS SUPPLIED

- Cuffs 2 x
- USB cable (A connector to Micro B connector) 2 x
- USB mains adapter 1 x
- Operating instructions 1 x

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any films, stickers or transport protection. **Never remove the rating plate and any possible warnings!**

Key numbers are shown as follows: (1)

Picture references are shown as follows:  
(Picture A)

## AT A GLANCE

(Picture A)

- 1 Hook-and-loop fastener
- 2 Cuff
- 3 Air chambers
- 4 Connecting socket
- 5  Switch heating function on / off (incl. control lamp)
- 6  Switch on / off + increase intensity level
- 7 Control lamp:
  - Green – intensity level I
  - Red – intensity level II
- 8 USB mains adapter
- 9 USB cable

## Dear Customer,

We are delighted that you have chosen the **MAXXMEE Air Compression Massage Cuffs**. Simply put them on, adjust them to your individual size with the hook-and-loop fastener and switch on. The integrated air cushions immediately begin to exert gentle pressure on the muscles from below, which stimulates the circulation and helps to regenerate your tired arms or legs. There are 2 intensity levels and a heating function, which you can adjust to your personal requirements.

If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website: [www.dspro.de/kundenservice](http://www.dspro.de/kundenservice)

We hope you have a lot of fun with your  
**MAXXMEE Air Compression Massage Cuffs.**

## Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

## SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Direct current

Alternating current

Protection class II

Protection class III

Only use indoors.

Energy efficiency class VI

Switching power supply

Safe electrical isolation

## SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

**DANGER**—high risk, results in serious injury or death

**WARNING**—medium risk, may result in serious injury or death

**CAUTION**—low risk, may result in minor or moderate injury

**NOTICE**—may result in risk of damage to material

## INTENDED USE

- The cuffs are intended for massaging and stimulating and heating your arm and leg muscles, and are used for relaxation.
- The cuffs are not designed for medical or therapeutic purposes.
- The device is for personal use only and is not intended for commercial applications.
- Use the product only for the specified purpose and only as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- The following people may use the device only following consultation with a doctor:
  - people who have no feeling in their arms or legs
  - people with circulation problems
  - women who are pregnant or have just given birth
  - people with a pacemaker
- **Only use the device if you have no physical illnesses or ailments!** If you are unsure whether the device may be used, consult a doctor.
- Do not use the device near to or on injuries, swelling, skin irritations, infected skin, warts, fungal infections, abnormal changes in the skin or body parts with no feeling.

## SAFETY NOTICES

 **WARNING:** Note all safety notices, instructions, illustrations and technical data provided with this device. Failure to follow the safety notices and instructions may result in electric shock, fire and / or injuries.

- Children younger than three years of age must not use this device as they are unable to react to overheating.
- The device must not be used by young children over three years of age, unless the regulation and / or control units have been preset by a parent or guardian or the child has been sufficiently briefed on how to operate the regulation and / or control units safely.
- This device may be used by children aged 8 and over and by people with reduced physical, sensory or mental abilities or a lack of experience and/or knowledge if they are supervised or have received instruction on how to use the device safely and have understood the dangers resulting from failure to comply with the relevant safety precautions.

- Children must not play with the device.
- Cleaning and user maintenance must not be carried out by children without supervision.
- Children younger than 8 years of age and animals must be kept away from the device and the connecting cable.
- The device has a hot surface. People who are insensitive to heat must be careful when using the device.
- The device may not be used by people who are vulnerable and are unable to react to overheating.
- The device is not designed for use in hospitals.
- Do not use the device if it is wet or damp. It must be completely dry before it is used.
- No safety pins or other sharp or pointed objects must be attached to the device or inserted into it.
- Do not switch on the device when it is folded up. It must always be spread out prior to use to prevent the build-up of heat.
- The device must always be disconnected from the mains power when left unattended and before cleaning it.
- Always allow the device to cool down completely before it is washed, folded up and / or stored away.
- To ensure that the device is not damaged, fold it up loosely, do not press it together firmly or bend it. Do not place any objects, animals or people on the folded device.
- Only use the included USB mains adapter or proper USB mains adapters with the same technical specifications as the cuffs.
- Do not make any modifications to the device. If the device is damaged, it must be replaced by the manufacturer or customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards.

**DANGER – Danger of Electric Shock**

- Charge and store the device only in closed rooms.
- Do not operate, charge or store the device in rooms with high humidity.
- Never immerse the control unit or the connecting cable in water or other liquids! Make sure that the control unit and the connecting cable cannot fall into water or become wet.
- If the control unit becomes wet or falls into water while it is connected to

the mains power, pull the mains plug out of the plug socket immediately. **Only after this** should you retrieve the control unit from the water. Never attempt to pull it out of the water while it is connected to the mains power!

- Never touch the device or the connecting cable with wet hands when the device is connected to the mains power.

#### **WARNING – Danger of Fire**

- Keep the device and accessories away from fire and other heat sources.
- Do not cover the control unit during operation or charging in order to prevent the device from catching fire.

#### **WARNING – Danger from Rechargeable Batteries**

- The control unit with a built-in battery must not be taken apart, thrown into a fire or exposed to high temperatures of over +40°C. The battery could ignite or explode.
- Make sure that the control unit cannot become wet. Do not use the control unit if liquid has got inside it. The battery in the control unit could ignite or explode.
- Do not cover the device during charging to prevent overheating.

#### **WARNING – Danger of Injury**

- Keep the packaging material away from children and animals. There is a danger of suffocation.
- If the battery leaks, avoid physical contact with the battery acid. Wear protective gloves if necessary. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately. There is a danger of burns from battery acid.

- Make sure that the USB cable is always kept out of the reach of small children and animals. There is a danger of strangulation.
- Make sure that the connected USB cable does not present a trip hazard.

#### **NOTICE – Risk of Damage to Material and Property**

- The device must be switched off when it is connected to or disconnected from the mains power.
- Make sure that the connecting cable is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- Pull the mains adapter out of the socket if an error occurs during charging or before a thunderstorm.
- When the USB mains adapter is pulled out of the plug socket, always pull on the USB mains adapter and never the USB cable.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- Protect the device from heat sources, flames, sub-zero temperatures, persistent moisture, wet conditions and impacts.
- Do not place anything on the cuffs so that they are not damaged.

## CHARGING UP THE BATTERY

### Please Note!

- Do not charge up the device via a computer because the computer might be damaged by the high charging current.
- Fully unwind the USB cable before connecting the USB mains adapter.
- Lay the connecting cable in such a way that it does not become squashed or kinked and does not come into contact with hot surfaces.

It is essential that you charge up the battery fully before you first use the device in order to achieve the maximum charging capacity. If the battery is fully discharged during use, it must be fully charged up again before it is next used.

1. Connect the USB cable (9) to the connecting socket (4) of the control unit for a cuff (2) and to the connecting socket of the USB mains adapter (8).
2. Only connect the USB mains adapter to a plug socket that is properly installed and matches the technical data of the mains adapter. The plug socket must also be readily accessible after connection so that the connection to the mains can quickly be isolated. The control lamp (7) in the control unit lights up red. The battery is charging.
3. When the control lamp lights up green, the battery is fully charged.
4. First disconnect the mains adapter from the power supply and then the USB cable from the mains adapter and the cuff.

## USE

### Please Note!

- A tingling sensation / itching during or after the massage is normal due to the circulation-boosting effect.
- Switch off the device immediately if you experience pain or feel unwell during use.
- To prevent overheating, do not use the device under a blanket or similar cover.
- Do not use the device for longer than 20 minutes without a break.

Start with the slowest level and gradually increase the level of intensity.

1. Place the cuffs (2) around your arms, feet or legs (Picture B) and fasten them with the hook-and-loop fasteners (1).  
Make sure that they are not applied too tightly so that they do not become too tight during operation.
2. Press (6) on the control units to activate the massage. The control lamp (7) lights up green. Intensity level 1 is switched on.  
The air chambers in the cuffs are slowly inflated in two cycles and emptied again after a few seconds.
3. If required, set intensity level 2 by pressing again. The control lamp lights up red.  
**Please Note:** When switched off and on again, the massage starts again at intensity level 1.
4. If required, also activate the heating function by pressing (5). The control lamp in lights up red.
5. After 10 minutes the massage stops automatically. It can be stopped at any time by pressing .

 The device should not be used before going to sleep. The massage has a stimulating effect, which may lead to problems getting to sleep.

## CLEANING AND STORAGE

### Please Note!

- Do not use any caustic or abrasive cleaning agents or cleaning pads to clean the device. They can damage the material.
- Wipe down the cuffs (2) with a damp, soft cloth and a little mild cleaning agent if required. Then dry them thoroughly.
- Store the device and accessories in a dry, frost-free place that is protected from direct sunlight and dust and cannot be accessed by children or animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

Do not attempt to repair an electrical device yourself!

**The cuffs cannot be switched on.  
The control lamp (7) flashes red.**

The battery is empty.

- ▶ Charge up the battery.

**No heat is produced.**

The device is not switched on.

- ▶ Press  (6). The control lamp (7) ought to light up.

The heating function is not switched on.

- ▶  Press (5). The control lamp in the switch ought to light up.

**The battery is not charging. The control lamp (7) in the control unit does not light up.**

The USB cable (9) is not plugged correctly in the connecting sockets on the control unit or the USB mains adapter is not plugged correctly in the plug socket.

- ▶ Correct the fit of the cable or adapter.

The plug socket is defective.

- ▶ Try another plug socket.

There is no mains voltage.

- ▶ Check the fuse of the mains connection.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of the cuffs with a built-in battery as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

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All rights reserved.
Gefabriceerd voor DS Produkte GmbH