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## LIST OF CONTENTS

Included in delivery	8
Product overview	8
Symbols	9
Signal Words	9
Intended use	9
Safety instructions	10
Use	12
Cleaning and Storage	12
Troubleshooting	13
Disposal	13
Technical Data	14

## INCLUDED IN DELIVERY

- 1 x massage device
- 1 x mains adapter
- 1 x operating instructions

Check the scope of delivery for completeness and the components for transport damage. If you find any damage, do not use the product but contact our customer service department.

Remove any possible films, stickers or transport protection from the product.

**Never remove the rating plate and any possible warnings!**

## PRODUCT OVERVIEW

(Picture A)

- 1 2x2 massage balls, heating and red light element
- 2 Control Panel
- 3 Hand loops
- 4 Connecting socket for the mains adapter
- 5  HT button: Red light and heat can be switched off or on again
- 6  On/Off button: When switched on the massage device starts the massage and changes direction approx. every minute. The red light and heat function is automatically switched on at the same time.
- 7  DR change of direction button: The massage balls massage continuously either clockwise or anticlockwise.

**Dear Customer,**

Thank you for choosing our **MAXXMEE** massage device.

If you have any questions about the device and about spare parts/ accessories, contact the customer service department via our website: [www.ds-group.de/kundenservice](http://www.ds-group.de/kundenservice)

We hope you have a lot of fun with your **MAXXMEE** massage device.

**Information About the Operating Instructions**

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

**SYMBOLS**

Danger symbols: These symbols indicate risks of injury. Read through the associated safety instructions carefully and follow them.

Supplementary information

Read operating instructions before use!

Circuit symbol for direct current

Protection class II

Protection class III

Energy efficiency class V

May only be used indoors!

Polarity of the coaxial connector (mains adapter)

**SIGNAL WORDS**

Signal words indicate hazards if the associated notices are not followed.

**DANGER** – high risk, results in serious injury or death

**WARNING** – medium risk, may have serious injury or death

**CAUTION** – low risk, may result in minor or moderate injury

**NOTICE** – may result in risk of damage to material

**INTENDED USE**

- This device is intended to be used for loosening the neck and back muscles and delivers relaxation. The device may also be used to massage abdominal, buttock and leg muscles. It is not a replacement for medically required treatment.
- The device is for personal use only and is not intended for commercial applications.
- The device should only be used as described in these instructions. Any other use is deemed to be improper.
- The following are excluded from the warranty: all defects caused by improper handling, damage or unauthorised attempts at repair. This also applies to normal wear and tear.

## SAFETY INSTRUCTIONS

**WARNING: Note all safety notices, instructions, illustrations and technical data provided with this device.**

Failure to follow the safety notices and instructions may result in electric shock, fire and/or injuries.

- This device should not be used by people (including children) with impaired sensory or mental faculties or who do not have adequate experience and/or knowledge of such devices, unless they are supervised by a person responsible for their safety, or have received instructions from such a person on how to use the device.
- Children must be supervised to ensure that they do not play with the device.
- The device has a heating function. People who are sensitive to heat must be careful when using the heating function.
- Keep sharp and pointed objects away from the device to ensure that it is not damaged.
- The device must always be disconnected from the mains power when it is left unattended and before cleaning it.
- In the interest of your own safety, you should inspect the device for damage each time before you use it. It may only be used if it has no visible signs of damage!
- Do not make any modifications to the device. Do not replace the connection yourself. If the device or the connection are damaged, they must be replaced by the manufacturer, customer service department or a similarly qualified person (e. g. specialist workshop) in order to avoid any hazards.

### Health-Specific Instructions

- Only use the device if you have no physical illnesses or ailments.
- The following people may use the device **only after consultation with their doctor**:
  - people who wear a pacemaker or similar medical implants;
  - people with heart disease or who suffer from epilepsy;
  - people who have back complaints or spinal disorders;
  - people who suffer from osteomyelitis, bone inflammation (in particular if the spinal column is damaged) or from osteoporosis;
  - women who are pregnant or have just given birth;

- people who suffer from circulation problems or thrombosis;
- people who have undergone an operation for which the treatment is not yet complete;
- people with physical abnormalities;
- people who have a malignant tumour.
- If you are unsure whether the device may be used, consult a responsible doctor.
- Do not use the device near to or on injuries, swelling, skin irritations, infected skin, warts, fungal infections, abnormal changes in the skin or body parts with no feeling.
- The device must not be used for longer than 30 minutes continuously. Do not massage one group of muscles for longer than 15 minutes.
- If you experience an unpleasant sensation while you are using the device, stop performing the massage immediately!
- Do not fall asleep while you are using the device!
- Do not use the device before going to bed. The massage has a stimulating effect, which may lead to problems getting to sleep.

### **DANGER – Danger of Electric Shock**

- Use and store the device only in closed rooms.
- Do not operate or store the device in rooms with high humidity.
- Never immerse the device or the connection in water or other liquids and ensure that they cannot fall into water or become wet.
- Should the device fall into water, switch off the power supply immediately. Do not attempt to pull the device out of the water while it is still connected to the mains power!

- Never touch the device or the connection with wet hands when the device is connected to the mains power.

### **WARNING – Danger of Injury**

- Keep plastic film and bags well away from children and animals. There is a danger of suffocation.
- Make sure that the connection is always kept out of the reach of small children and animals. There is a danger of strangulation.
- Do not insert any fingers or other body parts into the device. Do not open the zip during use.
- Make sure that the connected mains cable does not present a trip hazard. There is a danger of falling.

### **WARNING – Danger of Fire**

- Do not use the device under a blanket or similar item or in the vicinity of combustible material.

### **NOTICE – Risk of Damage to Material and Property**

- Only connect the mains adapter to a plug socket that is properly installed and matches the technical data of the device. The plug socket must also be readily accessible after connection so that the connection to the mains can quickly be isolated.
- Only use suitable extension cables whose technical data is the same as that of the device.
- In the event of damage / faults (e.g. unusual noises, odour or smoke development) and in a thunderstorm, switch the device off immediately and pull the mains adapter out of the plug socket.
- The device must be switched off when inserting or removing the

mains adapter. Always pull the mains adapter and never the mains cable to disconnect the device from the mains power.

- Do not lie or stand on the device.
- Protect the device from heat, naked flames, direct sunlight, extreme temperatures, moisture and collisions.
- Keep the device in a dry place at temperatures of between 0 and 40°C.

## USE

### Please Note!

- No pressure should be exerted on the device as this will slow down or stop the massage and the device could be damaged. Do not stand or lie on the device.
- Do not use the device for longer than 30 minutes without stopping.
- Do not massage one group of muscles for longer than 15 minutes.
- For reasons of hygiene, avoid direct contact of the skin with the device wherever possible. It should not come into contact with skin that has had cream or oil applied to it.
- A tingling sensation/itching during or after the massage is normal due to the circulation-boosting effect.

We recommend that you wait an hour after mealtimes before using the device.

### Areas of Application

When this device performs a massage, two massage heads massage the muscles twice. This is intended to relax them and boost circulation. The device is designed in particular for use on the upper back area, in particular the neck (Picture B). But it can also be used to

massage abdominal (Picture D), buttock (Picture E) and leg (Picture F,G) muscles.

### Operating the Device

1. Connect the mains adapter to the connecting socket of the device.
2. Insert the mains plug into a plug socket that is easily accessible.
3. Place the massage surface of the device onto the part of the body which is to be massaged. For a shoulder or back massage, place the straps of the device around your shoulders and fix them in place by putting your arms through the loops and thus holding the device in the correct position (Picture H).
4. Start the massage by pressing the  button.
5. After 15 minutes the device automatically switches off. To switch it off sooner, press the  button.
6. Pull out the mains plug and store the device in a dry place that is not accessible to children.

## CLEANING AND STORAGE

### Please Note!

- The device may not be connected to the mains and must be completely switched off for cleaning!
  - Do not use any caustic or abrasive cleaning agents for cleaning. These can damage the surface.
- Wipe down the device just with a damp cloth if necessary. Then wipe it down fully with a dry cloth.
  - Do not wind the mains cable around the device or the mains adapter.
  - Keep the device in its respective packaging in a dry place which is protected from sunlight and is not accessible to children and animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

### The device does not work.

The  button has not been pressed properly.

- ▶ Press the  button.

The mains adapter is not inserted correctly in the plug socket.

- ▶ Push the mains adapter into the plug socket as far as it will go.

The device is not connected correctly to the mains adapter.

- ▶ Correct the connection.

The plug socket is defective.

- ▶ Try another plug socket.

There is no mains voltage.

- ▶ Check the fuse of the electrical distributor (fuse box).

There is too much weight on the device.

- ▶ Reduce weight on the device.

### No heat is produced.

The red light and heat are not switched on.

- ▶ Press the HT button.

### The device gets too warm.

The device was in operation several times in succession for 15 minutes / too long in each case.

- ▶ Pull the mains adapter out of the plug socket and allow the device to cool down. Always disconnect the device from the mains when it is not in use.

### The device makes noises.

Operating noises are a normal consequence of the mechanics and not a defect.

### Smoke emerges from the device.

The device is defective.

- ▶ Pull the mains adapter out of the plug socket immediately and contact the customer service department.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

The adjacent symbol (crossed-out bin underlined) means that used devices should not be disposed of in the domestic waste, but should instead be taken to special collection and return systems.

Owners of used devices from private households can take them free of charge to the collection points of the public waste disposal authorities or the disposal points which are set up by manufacturers and distributors under the Electrical Equipment Act so that they can be disposed of in an environmentally friendly way and valuable raw materials can be recovered. If they are not disposed of properly, toxic ingredients may escape into the environment and cause harm to the health of humans, animals and plants. Businesses which market electrical and electronic devices are also obliged to take back used devices.

Waste prevention makes an even more valuable contribution to protecting the environment. This means if possible, as well as continuing to use or repair a device yourself, passing it on to another user also represents an ecologically beneficial alternative to disposal.

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All rights reserved.