Anleitung / Manual

KLIMMI SET

Art.Nr. 2873507

www.artzt.eu

ARTZT vitality Klimmi

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www.artzt.eu

ARTZT vitality Klimmi

is a practical tool for a connective tissue massage. Benefit: Adhesions and tensions of fasciae can be detected and treated with different intensity.

The wrists and finger joints are relieved and protected. Thanks to the soft silicone and low weight, you retain full control over the required pressure.

The targeted and punctual treatment with the Klimmi is an optimal supplement to the fascial roll.

Application areas

Klimmi was developed for the therapist's daily treatment of patients, but can also be used for self-treatment by the end user or patient.

It can be used to treat all areas of the body, especially the leg muscles (e.g. after recurring cramps and after sports), feet (e.g. during regular standing activities), forearms and ball of the thumb (e.g. during regular computer work).

General information

Apply sufficient oil to the skin before using the Klimmi. The Klimmi will thus glide evenly over the areas to be treated. The hand holding the Klimmi during treatment should remain dry to prevent the tool from slipping.

The flat side of the Klimmi always lies in the hand, the „nose“ builds up pressure in the tissue to be treated.

The movements during the treatment should be slow and even.

Slight pressure pain may occur during application. However, this should be perceived as quite bearable. If the pain is too severe, discontinue treatment or reduce the pressure.

ARTZT vitality Klimmi is not a toy!

Application

Klimmi can be used both lengthwise and crosswise. If you use Klimmi with the „nose“ in front, trigger points and small area adhesions can be relaxed. (Image 1)

If the Klimmi is held horizontally, it is suitable for the treatment of flat adhesions and muscle cramps. (Image 2)

1)

2)

When self-applying, the direction of treatment should be away from the centre of the body. That way, the arm's natural movements and thus strength are preserved.

In the case of mutual massage (e.g. in hard-to-reach areas such as the back), the direction of treatment can be varied as needed.

Degrees of hardness / colours

Soft / purple:
For small and sensitive areas such as the fore arm and nape of the neck

Hard / grey:
For large areas such as thighs and calves as well as deep-seated tissue structures

Material & care

ARTZT vitality Klimmi is made of silicone and is therefore absolutely skin-friendly and resistant to soap suds and oils.

If needed, the Klimmi can be easily washed with soap and water or cleaned with alcohol-free disinfectant.

ARTZT vitality

Videos & techniques

Further information about ARTZT vitality Klimmi available on: www.artzt.eu/klimmi

Videos about the application techniques and body areas can be found on our YouTube channel:

www.youtube.com/artztgmbh

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www.artzt.eu

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©2018 Ludwig Artzt GmbH. Alle Rechte geschützt.
©2018 Ludwig Artzt GmbH.