**KENDOX**

# MUSCLE RELIEF

Hochfrequenz Massagepistole

Gebrauchsanleitung.....	05
Instructions for Use.....	06
Mode d'emploi.....	08
Manuale di utilizzazione.....	10
Gebruiksaanwijzing.....	12
Használati utasítás.....	14
Návod k použití.....	16
Návod na použitie.....	17
Instrucțiuni de utilizare.....	19

#1

#2

#3

DE

#1

**MESSAGEAUFSÄTZE**

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#1

**MESSAGE ATTACHMENTS**

- A. Ball head**  
For massaging arms, waist, thighs and other muscle groups.
- B. Flat head**  
For relaxation of the soles of the feet and large muscle groups.
- C. Pointed attachment**  
For targeted application.
- D. Forked head**  
For areas around the neck, spine and Achilles tendon.

#2

**CORDLESS DEEP TISSUE MASSAGE DEVICE**

- A. Round standard attachment
- B. High-intensity vibration core transmits maximum vibration to the body
- C. Powerful motor with up to 3,900 rpm
- D. Cordless and lightweight
- E. 6 intensity speeds with LED display
- F. Ergonomic handle for secure control

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EN

## SAFETY NOTICES

**Thank you for purchasing this device. Please read the entire instruction manual carefully before use. Retain these instructions for future reference.**

- This device should only be used by children aged 8 years and over, as well as by persons with reduced physical, sensory or mental abilities or lack of experience and knowledge, if they are supervised or have been instructed on the safe use of the device and understand the dangers associated with its use. Children must never be allowed to play with the device. Cleaning and user maintenance must never be performed by children without supervision.
- Check the device and its parts thoroughly for any evidence of damage.
- Do not attempt to use the device if parts are damaged or if the device no longer works. Never attempt to repair the device yourself; contact a qualified service representative for repairs.
- Use the device as described in the manual. Any misuse or incorrect handling can lead to malfunctions.
- The device is not intended as a substitute for professional medical advice, diagnosis, or treatment. Use of this product is at your own risk.
- Please consult a doctor before use if your health is impaired by afflictions such as inflammatory/allergic reactions, during pregnancy or breastfeeding, etc.
- Do not use as a pain relief aid. Do not massage any parts of the body that are swollen, inflamed, or have a rash.
- Only ever use on your body when it is dry. The device can also be used over clothing: press lightly and move around the body for about one minute per body area.
- Do not use on the head or upon bony areas of the body
- Only use the device on the soft tissues of the body such as muscles; cease use immediately if it causes pain or discomfort.
- Keep hair, fingers and other parts of the body away from the shaft and the back of the applicator as they may be crushed or become trapped.
- Bruising can occur if strong pressure is applied. Frequently check the area being treated, and stop at any sign of discomfort.
- Do not use the device continuously for more than 15 minutes, otherwise the overheat protection may be activated, and the device will stop functioning automatically.
- Do not block the ventilation opening or insert any objects into them.
- Do not immerse the device in liquids or spray liquids into the ventilation openings.
- Use only the accessories supplied with the device.
- Keep the device away from liquids and heat sources.
- Do not leave the device unattended while it is charging or in operation.

### **Caution:**

If you have any concerns or uncertainties, consult your doctor before using this product. Read the instruction manual carefully. Only use this product for its intended purpose. Manufacturers and distributors are not liable for any misuse. This device is not intended as a substitute for professional medical advice, diagnosis, or treatment. Do not massage areas of the body that are swollen, inflamed, or have a rash.

### **CHARGING**

1. The average operating time of the battery is around 1.5 hours, but this

- is dependent on the speed set and the pressure applied during use.
2. Please charge the device fully before using it for the first time.
  3. Connect to the battery's charging port, and plug the adapter into the wall socket. When the red indicator light turns green, charging is complete.

#### THE DEVICE MUST NOT BE USED WHILE IT IS CHARGING.

4. Disconnect the device from the adapter after each charge.

#### OPERATING STEPS

1. Choose the appropriate detachable head, and press firmly into the massager device.
2. Press the power button, the red light will come on, and the device will start in Mode 1.
3. Press the power button again to set the optimal speed. You can select from Mode 1 to Mode 6.
4. Massage the desired body part with the most comfortable intensity, applying pressure as needed.
5. To switch off the device, press and hold the power button for 3 seconds.
6. If the device's performance decreases, you should charge it.

#### RECOMMENDATIONS FOR USE

- We recommend - depending on the desired effect - using the device for between 15 seconds 2 minutes per muscle group.
- A full-body treatment should not exceed 15 minutes (the device switches off automatically after 15 minutes of continuous operation to protect against overheating).
- The device can be used upon the same muscle group 2 to 3 times per day.
- Adjust the device to as gentle or as strong a setting as you feel comfortable with. If you feel that the massage is too weak, increase the intensity. If you feel any discomfort or pain as a result of application, select a lower intensity level.

#### RECOMMENDATIONS FOR USE

- We recommend - depending on the desired effect - using the device for between 15 seconds 2 minutes per muscle group.
- A full-body treatment should not exceed 15 minutes (the device switches off automatically after 15 minutes of continuous operation to protect against overheating).
- The device can be used upon the same muscle group 2 to 3 times per day.
- Adjust the device to as gentle or as strong a setting as you feel comfortable with. If you feel that the massage is too weak, increase the intensity. If you feel any discomfort or pain as a result of application, select a lower intensity level.

#### 1. Muscle activation

We recommend using the device for up to 30 seconds before training or physical activity, in order to achieve rapid activation of the muscles involved.

#### 2. Muscle reactivation

For reactivating the muscles during training or physical activity. Apply to the target muscle group for 15 seconds. This re-vitalizes the muscles and keeps them active.

#### 3. Recovery & regeneration

Use the device on the affected muscle group for a maximum of 2 minutes. We recommend application at 30-second intervals, with a 5 to 10 second break between applications.

#### 4. Relaxation

Use the device on a muscle group for 30 seconds to relax it. This is usually sufficient time to allow the muscle to relax.

#### MAINTENANCE

1. Clean with a soft dry cloth.
2. Do not use the device for more than 15 minutes at a time. It is recommended that you wait half an hour before starting to use again.
3. If the battery is not going to be used for a long time, charge and discharge it intermittently.

SPECIFICATIONS	
Battery capacity	2000 mAh
Power supply via DC adapter	110-240V/50-60Hz 0.5A
Output voltage	8.5V 1A
Power consumption	35W
Speed in Mode 1	2900 rpm
Speed in Mode 2	3100 rpm
Speed in Mode 3	3300 rpm
Speed in Mode 4	3500 rpm
Speed in Mode 5	3700 rpm
Speed in Mode 6	3900 rpm
Charge time	2.5 hours
Running time (after a full charge)	1.5 hours

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**CH:** MediaShop AG | Industriering 3 | 9491 Ruggell | Liechtenstein
**EU:** MediaShop GmbH | Schneiderstraße 1, Top 1 | 2620 Neunkirchen | Austria