 DE - Sicherheitshinweise	2
 EN - Safety Instructions	3
 IT - Istruzioni di Sicurezza.....	4
 FR - Consignes de Sécurité	6
 ES - Instrucciones de Seguridad.....	7
 SE - Säkerhetsanvisningar	8
 PL - Instrukcje Bezpieczeństwa	10
 NL - Veiligheidsinstructies.....	11

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EN - Safety Instructions

Declaration of Conformity

The ION-Sei electric sonic toothbrush complies with the requirements of the European Directives 2014/30/EU (Electromagnetic Compatibility) and 2011/65/EU (Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment).

Warnings

- Although the electric sonic toothbrush is not inherently dangerous, you should be aware of the potential risks associated with its use and act in a way that prevents harm to yourself and others. Therefore, please follow these rules:
- Only operate the electric sonic toothbrush in environments that meet the specifications in the Technical Data section on page 14. High temperatures and humidity can damage the built-in battery and other components, potentially causing a fire or explosion. The mains voltage must match the information on the charging device's nameplate.
- The charger is not waterproof. Do not use the charger outdoors. Do not touch the charger with wet hands. Ensure that the charger is always dry and does not come into contact with water.
- The electric sonic toothbrush is not a toy. Children over the age of 10, as well as persons with reduced physical, sensory, or mental abilities, or lack of experience and knowledge, may use the electric sonic toothbrush only if supervised or instructed on its safe use and are aware of the potential risks.
- Do not use the electric sonic toothbrush if you are a pacemaker user or have any other implant.
- If you are undergoing dental treatment, consult your dentist before using the electric sonic toothbrush. If you experience excessive gum bleeding after use, consult your dentist.
- Unplug the charger before cleaning it and ensure it is completely dry before reconnecting it. Never immerse the charging station or the charger completely in water or any other liquid.

Disposal / Recycling

Please only pass on the ION-Sei electric sonic toothbrush along with the user manual. The toothbrush and charger contain high-quality materials and components that can be recycled and reused.

The symbol indicates that electrical or electronic equipment must be disposed of separately from normal household waste. In the European Union, there are various collection systems for used electrical and electronic products. Dispose of these marked devices at your local collection and recycling point.

Before disposing of the handle, remove the battery and dispose of it separately. The battery must not be disposed of with household waste.

If unsure or lacking the necessary tools, have the battery removed by your dealer or a professional.

⚠ CAUTION: This process cannot be reversed. Once the cap at the bottom of the handle is removed, the toothbrush should not be used.

- ▶ Remove the toothbrush from the charging base, turn it on, and let it run until the battery is completely drained.
- ▶ Unscrew the screw at the bottom of the handle using a suitable Phillips screwdriver and remove the cap.
- ▶ Detach the attached wire spiral and bend it to the side.
- ▶ Insert a small screwdriver into the gaps and gently rotate to unlock the two locks between the inner part and the housing.
- ▶ Hold the handle upside down and press the titanium dioxide rod vertically onto a table surface to push the inner part out of the housing towards the open side.
- ▶ Remove the battery and unplug the cable connector.

Intended Use

The ION-Sei electric sonic toothbrush is intended exclusively for cleaning and caring for human teeth. Any other use is expressly prohibited.

The ION-Sei electric sonic toothbrush is intended for private use only and is not for use by multiple patients in a dental practice or clinic. Using the device or parts of it in any way other than as described in this manual is considered misuse.

Further Information

If you have any questions about the use of the ION-Sei electric sonic toothbrush that are not answered in this manual, please contact your dealer.

For the latest information about the ION-Sei electric sonic toothbrush, visit www.ion-sei.com.

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