-  *Instruction Manual*
-  *Manual de instrucciones*
-  *Mode d'emploi*
-  *Bedienungsanleitung*
-  *Istruzioni per l'uso*
-  *Manual de instruções*
-  *Gebruiksaanwijzing*

International patents pending. All imitations will be prosecuted.
GYMFORM® ELECTRO FITNESS TRAINER is a registered EU/CTM trade mark.

SPORT

GYM[®]
form

ELECTRO FITNESS TRAINER


1

2


3




4

5

6

9

7

8

GYMFORM® ELECTRO FITNESS TRAINER INSTRUCTION MANUAL

Congratulations on your purchase of Gymform® ELECTRO FITNESS TRAINER. Your safety is our primary concern, so please ensure that you familiarize yourself with all the information in this manual before using your Gymform® ELECTRO FITNESS TRAINER. It is essential that anyone who intends to use this device reads and understands the instructions provided. An informed user is a safe user. Congratulations on your purchase of Gymform® ELECTRO FITNESS TRAINER. We hope you enjoy your sessions.

THE SET INCLUDES:

- 1 Gymform® ELECTRO FITNESS TRAINER unit with 6 programmes
- 1 Adjustable Gymform® ELECTRO FITNESS TRAINER vest
- 1 Adjustable Gymform® ELECTRO FITNESS TRAINER pants
- 2 Adjustable arm straps
- 20 Electrodes
- 1 Instruction manual
- 1 Sports bag
- 1 Diet booklet
- 1 Exercise guide
- 1 DVD Personal trainer guide

• WARNING

Important: After finishing a workout session, allow 48 hours to pass before doing another one. Before embarking on any exercise programme, including the Gymform® ELECTRO FITNESS TRAINER, consult your general practitioner and your physiotherapist.

Ensure that you read the instructions, warnings and precautions carefully before using your Gymform® ELECTRO FITNESS TRAINER.

Gymform® ELECTRO FITNESS TRAINER is designed for healthy adults. However, the device should be used with care. Read the following guidelines and instruction manual before using the product.

- Using the Gymform® ELECTRO FITNESS TRAINER while connected to high-frequency surgical equipment may cause burns at the contact points of the Gymform® ELECTRO FITNESS TRAINER and may result in damage to your Gymform® ELECTRO FITNESS TRAINER.
- Operation of a short-wave therapeutic device or microwave oven in the general vicinity (e.g. 1 metre away) may cause instability in the output power of the Gymform® ELECTRO FITNESS TRAINER.
- This is a powerful electrical muscle stimulation device and may cause injury if you start with the most advanced exercise programme right away. Start with a lower level and gradually increase the intensity.
- Do not use on animals.
- Do not use while sleeping.
- If you observe damage or imperfections in the device or it gives off a strange odour, makes noise or produces smoke, do not turn the device on.
- Do not attempt to open, tamper with or repair the device, as that will render the warranty null and void.
- Use only Gymform® ELECTRO FITNESS TRAINER components.
- Do not use in the bathroom or in other places with a high amount of moisture or humidity.
- The batteries should be removed from the device before it is disposed of.
- The batteries should be safely disposed of.
- Stop the device immediately if you feel ill, dizzy or pain in the chest.

Do NOT use Gymform® ELECTRO FITNESS TRAINER:

- If you have an electronic implant (such as a pacemaker) or suffer from another heart problem.
- If you are, or think you might be, pregnant.
- If you have cancer, epilepsy or are under medical supervision for a cognitive disorder.
- If you are breastfeeding.
- If you have a high fever or an acute bacterial or viral disease.
- If you have a femoral or abdominal hernia.

- If you have a disease of the nervous system.
- If you have a dermatological disease.
- If you have multiple sclerosis.
- If you have symptoms of acute or chronic disorders in general.

Warnings for women: Before using Gymform® ELECTRO FITNESS TRAINER, you should wait:

- At least six months after the birth of your child (after consulting your doctor).
- A month after having an intrauterine device (IUD) fitted.
- At least three months after having a caesarean section (after consulting your doctor).
- Until after the first days of your menstrual period.

Before using Gymform® ELECTRO FITNESS TRAINER, you should obtain permission from your doctor or physiotherapist:

- If you suffer from back problems. In this case, the device should only be used at low intensity levels.
- If you suffer from a serious disease or injury not mentioned in this section.
- If you were operated on recently.
- If you are diabetic and take insulin.

When placing the electrodes, always remember:

- Place the electrodes only on the muscles groups they are designed for. Ensure the entire surface of the electrode is in contact with the skin. Otherwise the current density may be increased, which could cause pain and damage the skin.
- Do not connect, disconnect or tamper with the electrodes during the workout. First turn off the device.

- Do not place electrodes on recent scars, inflamed skin, burns, ulcers, infected areas or areas that are susceptible to acne, thrombosis or other vascular problems (e.g. varicose veins), moles or any part of your body with limited sensitivity.
- Do not place on areas with injuries or limited mobility (e.g. fractures or sprains).
- Do not place directly on metal implants.
- Wipe after each use.

Possible Adverse Reactions:

There have been isolated reports of skin reactions caused by the use of electro-stimulation devices, including allergies, prolonged skin redness and acne.

If any of these symptoms occur, switch off the device and consult your doctor.

If you have never exercised before or have not exercised for a long time:

Consult your general practitioner.

Precautions:

- Do not use Gymform® ELECTRO FITNESS TRAINER if you are drowsy or while sleeping.
- Do not use Gymform® ELECTRO FITNESS TRAINER if your body is in contact with water.
- Do not use Gymform® ELECTRO FITNESS TRAINER in combination with any other device that transmits electrical current to your body or with any other machine.
- Do not connect the metal snap fasteners or electrodes to any other object.
- Use only the electrodes supplied by Industex.
- Do not attach two electrodes together.
- If any of the electrodes are not attached properly, you will notice an itching sensation.
- Do not wash the suit in the washing machine. Follow the washing instructions (see below).
- If suitable supervision and instructions are provided, the device may be used by children aged 8 and older, by people with reduced physical, sensory or mental abilities, and by people lacking in experience and knowledge.

- Children should not be allowed to play with the device.
- The cleaning and maintenance to be performed by the user should not be done by children without supervision.

PRODUCT DESCRIPTION:

1. Gymform® ELECTRO FITNESS TRAINER unit
2. Large rectangular electrodes
3. Small rectangular electrodes
4. Round electrodes
5. Vest
 - 5.a) Adjustable armbands
 - 5.b) Adjustable top bands
 - 5.c) Vest cable
 - 5.d) Adjustable torso straps
6. Pants
 - 6.a) Legs
 - 6.b) Pant belt
 - 6.c) Strips with quick opening and closing system for the pants
 - 6.d) Pants cable
7. Console connection for pants
8. Display
9. ON/OFF button 
10. UP button
11. DOWN button
12. MASTER button
13. SET button
14. STOP button
15. Programme
16. Intensity of different muscle groups
17. Workout time remaining
18. Battery symbol
19. Console connection for vest
20. Snap
21. Reset button

Before performing the following steps, make sure the unit (1) is in the OFF position.

Battery Installation:

- Open the battery compartment cover (batteries not included) on the back of the unit (1) and insert four LR20 alkaline batteries.
- Insert the four batteries in accordance with the polarity marks indicated on the battery compartment cover.
- Replace the battery compartment cover.

Electrode placement (see Figure 1): Use the snaps (20) to attach the electrodes to the suit. Each electrode has one or two snaps on the back. To attach the electrode to the suit, snap it into place until it clicks (See Figure 4). There are three kinds of electrodes and they should be placed as follows:

- The four small rectangular electrodes (3) should be placed in the abdominal area of the vest (5) (stomach area).
- The four round electrodes (4) should be placed on the adjustable armbands (5.a).
- The large rectangular electrodes (2) should be attached to the snaps in the area of the chest, the trapezius, the lumbar area, the glutes, and the front and back of the thighs (see Figure 1).

Electrode conductivity:

To ensure that Gymform® Electro Fitness Trainer works properly, it is necessary to ensure electrode conductivity. Choose the method that works best for you from the following:

- Apply Gymform® conductive gel on each electrode. Gymform® gel is available at the same establishment where you purchased the Gymform® Electro Fitness Trainer®.
- Use a wet sponge or cloth to moisten the part of each electrode that makes contact with the skin.
- You can also wear a damp T-shirt and pair of shorts under the Gymform® Electro Fitness Trainer suit. We recommend wearing garments made of cotton.

Putting on the Vest (see Figure 6):

Once the electrodes are in place and their conductivity has been ensured (see explanation in the Electrode Placement section), put on the vest (5). Close the zip. The suit has to fit snugly to work well. First pull on the two adjustable top straps (5.b). Then continue moving down and adjust the adjustable torso straps (5.d). Pull on them to adjust the vest, as shown in Figure 6, and make sure that all the electrodes are in contact with the skin. The vest (5) should fit snugly, but comfortably. Then put on the armbands (5.a). Make sure the armband (5.a) is firmly snapped to the vest (5). Open the strip with the quick opening and closing system on the adjustable armband. Place the armband around the upper arm (biceps area) and attach it by threading the strap through the plastic strap adjuster and attaching it with the quick opening and closing system strip (See Figure 8). Ensure it is firmly attached and the electrodes are in contact with the skin. The proper placement of the armbands (5.a) is shown in Figure 9.

Putting on the Pants:

1. While the pants (6) are completely open, hold them by the pant belt (6.b) and wrap the pants (6) around the back of your waist. The place for electrodes on the back should be in the area of the glutes. Attach the pant belt (6.b) with the quick opening and closing system strip.
2. Then grab one pant leg (6.a), wrap it around your thigh and attach it with the quick opening and closing system strips for pants (6.c). The pant leg (6.a) should be wrapped around the very top part of the thigh. If it is too high, the problem is that the pant belt (6.b) is not in the right place

and should be lowered.

- Do the same thing with the other pant leg (6.a).
- When done correctly, there are two electrodes on the glutes, two electrodes on the front part of the thighs and two electrodes on the back part of the thighs. Now make a final check to adjust all the quick opening and closing system strips and ensure that the pants (6) are perfectly adjusted to your skin. You should feel comfortable.

Connecting the Suit to the Gymform® ELECTRO FITNESS TRAINER unit (see Figure 2):

- Make sure the unit (1) is turned off.
- Connect the pant cable (6.d) into the unit (1) via the console connection for pants (7).
- Connect the vest cable (5.d) into the unit (1) via the console connection for the vest (19).

Using the Gymform® ELECTRO FITNESS TRAINER unit (see Figure 3)

The unit menu is intuitive. The more you use the device, the easier it will be.

- Press the ON/OFF button (9) to turn on the unit (1). The display (8) will come on.
- The Gymform® ELECTRO FITNESS TRAINER unit (1) has the following functions to choose from:
 - PROGRAMME MODES (15): There are six different programme modes. The word MODE is shown on the display (8).
 - INTENSITY LEVELS (16) of eight different muscle groups.
 - WORKOUT TIME (17).

In the lower left of the display (8), there is a battery symbol (18) that indicates the remaining battery life. When the symbol is flashing, you should change the batteries very soon because they are about to go dead.

3. BUTTON FUNCTIONS:

- ON/OFF button (9): Used to turn the unit on and off.
- UP BUTTON (10): Used to move up on the display (8) and to increase the amounts in specific functions (Programme, Intensity Level and Time)
- DOWN BUTTON (11): Used to move down on the display (8) and to decrease the amount in specific functions (Programme, Intensity Level and Time)
- MASTER BUTTON (12): Used to change the intensity of all muscle groups at the same time.
- SET BUTTON (13): Used to select and validate each function chosen on the display (8). Each time you choose a function, programme or level, press the set button to validate it and go on to the next option.
- STOP BUTTON (14): When this button is pressed, the device will stop emitting current. If you feel uncomfortable and want to stop the exercise immediately, press this button. All the intensity levels will go back to 0.

PROGRAMME SELECTION

- Use the UP button (10) and the DOWN button (11) to move up and down the display (8) until you select the MODE function.
- Then press the SET button (13) to select this function.
- You'll notice that the numbers to the right of the word MODE will start to blink. Use the UP button (10) and the DOWN button (11) to choose the desired programme.
- Then validate the chosen program by pressing the SET button (13).
- You can choose from the following programme modes:

Programme 1	General Conditioning
Programme 2	Abdominals
Programme 3	Cardiovascular Exercise
Programme 4	Muscle Development
Programme 5	Massage and Drainage
Programme 6	Glutes and Thighs

Programme 1: General Conditioning:

Programme to become familiar with electrical muscle stimulation. This is a good level to start using the product. It is also designed for general body conditioning and maintenance. It will keep you in shape. Programme designed for use with the DVD Personal Trainer Guide (20-minute workout).

Programme 2: Abdominals:

This programme is designed to tone all the muscle groups, especially the abs.

Programme designed for use with the DVD Personal Trainer Guide (20-minute workout).

Programme 3: Cardiovascular Exercise:

Programmed designed to help slim down your figure.

Programme designed for use with the DVD Personal Trainer Guide (20-minute workout).

Programme 4: Muscle Development:

Programme designed to strengthen and increase your muscles. Programme designed for use with the DVD Personal Trainer Guide (20-minute workout).

Programme 5: Massage and Drainage:

Programme designed to use massage or drainage to relax your muscles after a workout session. Should be used 5 or 10 minutes after using a different programme.

Programme 6: Glutes and Thighs:

Programme designed to work on the glute and thigh area, as well as all the areas where fat accumulates, such as the stomach (ab area), the back and the chest. Use this programme with the exercises described in the Exercise Guide booklet.

Tip for beginners: The first few times you use the product, we recommend that you choose a program that provides a pleasant experience and a good sensation and that you do not pursue other benefits. One such option is Program 5: Massage

and Drainage. Start with low intensity levels at your own pace and follow the exercise instructions on the DVD or the Exercise Guide booklet. Become familiar with the device and gradually increase the intensity level if you wish. We recommend doing it in the first few sessions.

SELECTING THE INTENSITY LEVEL

After selecting the programme you want to work on, you have to choose the intensity level of the eight different muscle groups the Gymform® ELECTRO FITNESS TRAINER works on. As shown on the display (8) of the unit (1), the muscle groups are as follows (see Figure 5):

- Trapezius
- Chest*
- Arms
- Abdominals
- Lumbar
- Glutes
- Front of thigh
- Back of thigh

- Use the UP button (10) and the DOWN button (11) to move up and down the display (8) until you select the muscle group you want to work on.
- Select by pressing the SET button (13).
- When the number starts blinking, choose the desired intensity level by going up with the UP button (10) and going down with the DOWN button (11). As you increase the intensity level, you will feel the intensity in the different targeted muscle groups.
- When you have reached the desired intensity level, validate it by pressing the SET button (13). As you increase the intensity level, you will feel the intensity in the different targeted muscle groups.
- Repeat the process for all the muscle groups you want to work on.
- MASTER FUNCTION: When the MASTER button (12) is pressed, the levels for all the muscle groups blink at the same time and the intensity level of all of them can be raised or lowered at the same time. Tip: Once you have become familiar with the product. A practical approach is to start with the MASTER function and increase the intensity level of all the muscle groups at the same time. Then you can individually adjust the intensity depending on the sensitivity of each muscle group and the performance level you are seeking in each area.
- Now you can choose the workout time. The default workout time is 30 minutes, but you can set the number of minutes you want. Just 15 to 20 minutes a day is enough to achieve excellent results.
- Go to the TIME function by pressing the UP button (10) and the DOWN button (11). When the numbers start blinking, select using the SET button (13).
- Then choose the workout time using the UP button (17) and the DOWN button (11). Validate your choice by pressing the SET button (13).

* Women should pay special attention if they want to use the electrodes in the chest area, as the sensation can be unpleasant in this sensitive area. Start with very low intensity levels and stop if there is any feeling of discomfort.

Once you have chosen all the functions, you can start working out. You can choose to follow the exercises in the DVD Personal Trainer Guide (20-minute workout) or you can do the exercises explained in the Exercise Guide booklet.

•DVD Personal Trainer Guide (20-minute workout): Insert the DVD in the DVD player. Before starting the DVD, choose from programmes 1, 2, 3 and 4 on your Gymform® Electro Fitness Trainer unit. Then choose the desired intensity levels and the 20-minute workout time. Once the unit (1) is operating, start playing the DVD and do the same exercises in real time as the personal trainer is doing in the DVD. While you do the exercises, the Gymform® Electro Fitness Trainer suit will be working your muscles simultaneously to obtain the benefits of the chosen programme.

If at any time you find it is difficult for you to do an exercise on the DVD due to the high intensity level chosen for a specific muscle group, lower the intensity level to continue with the workout.

•Exercise Guide booklet: Read the booklet and choose the exercises you want to do. Follow all the instructions on each exercise.

During the workout session, you can modify any of the functions on the unit without stopping your workout. You can modify the intensity level, the workout time and the programme while the unit is still operating.

When changing programmes, make sure you lower the intensity levels first. Remember that a fast way to do it is to use the MASTER button (12). Otherwise, the intensities you have chosen for the current programme may not be appropriate for the new programme you choose.

Note:

Normal product performance may be altered by electrostatic disturbances. To return to normal operation, press the RESET button (21), which is inside the unit (1), by inserting a thin object into the hole and pushing.

AUTO OFF FUNCTION

If there is no activity for one minute, the display (8) will automatically shut off.

MEMORY FUNCTION

Each time you turn on the Gymform® ELECTRO FITNESS TRAINER unit, the display (8) will show the settings of the last

workout session, i.e. the same programme, same intensity levels of the muscle groups and the same workout time.

This function is very practical because it means you do not have to change the settings every time you want to work out. If you want to repeat the same workout session, press the SET button (13).

WARNING: If you decide to repeat the programme from the previous workout session, when you press the SET button (13), the programme will begin with the same settings as the last workout session. Gymform® ELECTRO FITNESS TRAINER will not supply the full intensity all at once, but will gradually increase the intensity in a few seconds' time. However, if the intensity level selected is too high, we recommend you start programming the workout from the beginning.

END OF WORKOUT

We recommend drinking a lot of water at the end of the workout. Besides hydrating the body, it helps to naturally eliminate burned-off fat with the help of the lymphatic system.

When you finish the workout session, turn off the unit. Disconnect the pants (6) and the vest (5) from the unit.

Loosen the adjustable straps as shown in Figure 7 and take off the vest (5) and pants (6).

We recommend cleaning the suit when you finish a workout session. See the Cleaning and Maintenance section below.

TECHNICAL SPECIFICATIONS:

Operating Voltage: 6.0V (SIZE D Battery x 4)

Operating Current: $\leq 2.1A$

Stand By Current: $\leq 0.1mA$

Operating Mode: 6 types

Output: 8 pcs

Output Voltage: Level 1: 2V VPP (Load 4K OHM)

Level 99: 100V VPP (Load 4K OHM)

Level 2~ Level 99 Every Type increase 1V VPP

This symbol indicates that the device is protected against electrostatic discharge by BF-type double insulation.

ISO 7000-1641 - User's guide

CLEANING AND MAINTENANCE

DO NOT WASH the vest (5) or the pants (6) of the Gymform® ELECTRO FITNESS TRAINER in the washing machine or by hand: JUST WIPE WITH A DAMP CLOTH AND ALLOW TO AIR DRY. The unit (1) can be cleaned with a damp cloth. Do not use products containing contaminants or acids, as they may damage the device.

Do not use products containing chemical abrasives to clean the electrodes.

To extend the life of the suit, follow the washing and maintenance instructions and use the suit with care. After each workout, the suit should be cleaned and dried and then stored with care. Store the suit with all the quick opening and closing system strips closed for the best maintenance.

PRECAUTIONS WITH BATTERIES

Ensure the batteries are properly installed before turning the unit on.

Turn off the unit to change the batteries.

Remove the batteries if you are not going to use the device for an extended period of time.

Take the necessary precautions when disposing of the batteries. Do not expose them to fire.

Remove the batteries from your Gymform® ELECTRO FITNESS TRAINER if you do not plan to use it for an extended period of time.

DISPOSAL OF MATERIALS

A symbol of a crossed-out wheeled bin means you should find out about and follow local regulations about disposing of this kind of product.

Do not dispose of this product as you would other household waste. Dispose of this device in accordance with the corresponding local regulations.

Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and/or human health and should be recycled properly.

The battery should be disposed of separately. To protect the environment, all batteries, including rechargeable ones, should not be disposed of as you would other household waste. Take them to a disposal centre and ensure it complies with applicable legislation on battery disposal.

Quality Warranty:

This product is covered by a warranty against manufacturing defects subject to the time periods stipulated by the legislation in force in each country.

This warranty does not cover damages resulting from inadequate use, negligent commercial use, abnormal wear and tear, accidents or improper handling.

Australia & New Zealand only

This product has the benefit of certain Consumer Guarantees. These are prescribed by the Australian Consumer Law & the New Zealand Consumer Guarantees Act 1993 both of which provide protection for consumers. There is no express warranty for this product in Australia or New Zealand. The above paragraph refers to other countries.

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ISL GmbH
Industex GmbH
Best Direct B.V.
Best Direct International Ltd,
ISL ITALY S.r.l.