Personenwaage Felicitas

Felicitas Bathroom Scale | Báscula personal Felicitas

Pèse-personne Felicitas | Bilancia pesapersona Felicitas

Waga osobowa Felicitas

Deutsch

English

Español

Français

Italiano

Polski

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BM 702 | 702-1 | 702-2 | 702-3 | 702-4

Operating instructions | Instrucciones de operación | Mode d'emploi

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Dear Customer!

You have chosen to purchase a high-quality product from the brand ADE, which combines intelligent functions with an exceptional design. These bathroom scales make it possible for you to determine and monitor your body weight. The many years of experience of the brand ADE ensures a high technical standard and proven quality.

We hope you enjoy weighing yourself and wish you success in reaching your goal weight!

Your ADE Team

General Information

About this manual

This operating manual describes how to safely use and care for your products. Keep this operating manual in a safe place in case you would like to refer to it later. If you pass this item on to someone else, pass on this manual as well.

Failure to observe the recommendations in this operating manual can lead to injuries or damage to your product.

Explanation of symbols

This symbol, combined with the word DANGER, warns against serious injury.

This symbol, combined with the word WARNING, warns against moderate and mild injuries.

This symbol, combined with the word NOTE, warns against material damage.

This symbol signifies additional information and general advice.

Table of Contents

General Information	8
Table of Contents	9
Intended Use.....	10
Safety	10
Overview (Scope of Delivery)	11
Calibrating	12
Weighing	12
Cleaning.....	12
Fault/Remedy.....	12
Technical Data	13
Guarantee	13
Disposal	13

Intended Use

The bathroom scale is designed for the weighing of people in the home. It is not suitable for use in clinics, medical practices, etc.

The maximum carrying capacity of 160 kg must not be exceeded.

Safety

- These scale can be used by children 8 years of age and above, as well as by persons with reduced physical, sensory, or mental capacities or lacking experience and knowledge, if they are supervised or have been instructed in the safe use of the scale and understand the dangers that may arise.
- Children must not play with the scale.
- Cleaning and user maintenance must not be carried out by children, unless they are 8 years of age or older and are being supervised.

DANGER for children

- Keep children away from packaging material. Risk of choking if swallowed. Also beware that children do not pull the packaging material over their heads.

WARNING against risk of injury

- Do not step onto the outer edge of any one side of the scale. Risk of tipping!
- Do not step onto the scale with wet feet or when its surface is damp. Risk of slipping!

BEWARE of material damage

- Place the scale on a stable, even surface.
- Do not overload the scale. Pay attention to the "Technical Data". Place yourself on the scale exclusively for weighing purposes.
- Protect the scale from sharp jolts and vibrations.
- Do not make any modifications to the scale. Only allow repairs to be carried out at a specialist workshop. Repairs carried out incorrectly can considerably endanger the user.

- There are feet on the underside of the scale. Since floor surfaces are made from various materials and are treated with diverse care products, the possibility cannot be completely excluded that some of these substances contain constituents that attack and soften the feet of the device. If necessary, place an anti-slip underlay beneath the scale.

Overview (Scope of Delivery)

- 1 Dial
- 2 Weight pointer
- 3 Weighing surface
- 4 Calibration wheel

Not depicted: Operating instructions

Calibrating

Before weighing yourself for the first time, you have to calibrate the scale in order to obtain a weighing result that is as exact as possible.

- Using the calibration wheel, centre the weight pointer on the "0" of the dial.

Check the zero-setting from time to time, e.g. when you change the location of the scale. When no load is on the scale, the weight pointer must be centred on "0".

Weighing

1. Position yourself on the scale. Remain standing still and do not hold on to anything.
2. Read off your measured body weight from the weight pointer.

- Weigh yourself with as little clothing on as possible.
- Try to weigh yourself always at the same time and under consistent conditions in order to obtain comparable values.

Cleaning

NOTE on material damage

- Do not immerse the scale in water or other liquids.
 - Do not under any circumstances use harsh, scraping, or abrasive cleaning agents or the like to clean the scale. If you do, the surface could be scratched.
- Wipe down the scale as needed using a lightly dampened, soft cloth.
 - Disinfect the weighing surface from time to time using a suitable disinfecting agent (available at specialist retailers and pharmacies).

Fault/Remedy

Fault	Possible Cause
The scale is showing obviously false weighing results.	<ul style="list-style-type: none">▪ Is the scale standing on a firm, level surface?▪ Is the scale correctly calibrated? (see chapter "Calibrating").

Technical Data

Model:	Felicitas Bathroom Scale BM 702 702: black, 702-1: mint, 702-2: white/grey, 702-3: red, 702-4: petrol
Measuring range:	max. 160 kg
Indicator division:	1 kg
Weight:	approx. 3 kg
Dimensions of scale:	approx. 275 x 345 x 145 mm
Distributor: (no service address)	Waagen-Schmitt GmbH Neuer Höltigbaum 15 22143 Hamburg, Germany

Our products are constantly being developed and improved.
For this reason, design and technical modifications are possible at any time.

Guarantee

The company Waagen-Schmitt GmbH guarantees for 2 years from the purchase date the cost-free rectification of defects on the basis of material or fabrication errors through repair or exchange.

In the event of a guarantee claim, please return the product with the proof of purchase (stipulating the reason for claim) to your dealer.

Disposal

Disposing of the packaging

Dispose of the packaging with similar materials. Put cardboard with waste paper, foil with recyclable materials.

Disposing of the product

Dispose of the product in accordance with the regulations that apply in your country.

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