**GYM**  
*form*®

**PERCUSSION MASSAGER**

-  *Instruction Manual*
-  *Manual de instrucciones*
-  *Mode d'emploi*
-  *Bedienungsanleitung*
-  *Istruzioni per l'uso*
-  *Manual de instruções*
-  *Gebruiksaanwijzing*

International patents pending. All imitations will be prosecuted. Gymform® Percussion Massager is a registered EU/CTM trade mark.

Gymform® Percussion Massager  
Percussion massager with rechargeable battery  
and four interchangeable heads.

## WARNINGS:

- Precautions: Read all of the instructions carefully before using the device for the first time.
- the appliance is only to be used with the supply unit provided
- This appliance contains batteries that are non-replaceable
- **WARNING:** For the purposes of recharging the battery, only use the detachable supply unit provided with this appliance
- this appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved
- children shall not play with the appliance
- cleaning and user maintenance shall not be made by children without supervision
- Do not immerse in water.
- Keep away from liquids and heat sources.
- Keep the air vent free of dust and debris.
- Do not remove any screws or disassemble the device.
- Do not use continuously for more than 30 minutes.
- Unplug the cable after charging or before using the device.

- Adult use only. Keep out of reach of children.
- NEVER leave the device unattended, especially if children are present.
- This device should NEVER be used by individuals suffering from a condition that may limit their capacity to operate it.
- Ensure that children do not use the device without adult supervision.
- NEVER cover the device while it is in use.
- NEVER use this device in bed.
- Do not use if you are injured. Ask your doctor before using this device.
- Do not use this device to complement or replace any medical treatment. Consult your doctor before using this device if you feel unexplained pain, are undergoing medical treatment or use any kind of medical device.
- This product is not for medical use. It is designed solely to provide a good massage.
- Consult your doctor before using this device if: You are pregnant. You have a pacemaker. You have any questions about your health.
- Not for use by diabetics or people with sensory deficiencies.
- NEVER use this device directly on rashes or swollen and inflamed skin.
- DO NOT use this device before going to sleep. Massage has a stimulating effect and may make it harder to fall asleep.
- Do not use on your head or near your genitals.
- Do not use:
  - . If you have an aneurysm, are bleeding, or take anticoagulants.

. If you have heart disease and wear a pacemaker or defibrillator.

. If you have cancer.

. Within 90 days following intra-articular fixation. The device should not be used within 10 cm of the site of internal fixation.

. On sensitive areas of the head, cervical vertebrae or spinal column.

. On areas near bones with less muscle, e.g. tibia, instep, back of the hand.

- It also comes with four heads with different functions:

o U-shaped head (2):

. For use on muscles near bones such as the spinal column.

. For use on muscles near tendons such as the Achilles tendon.

o Finger-shaped head (3): for targeting very specific areas with pain and tight muscles.

o Ball-shaped head (4): made of a soft material for use on large muscle groups and by inexperienced users.

o Flat head (5): covers larger surface areas. Can be used on all muscle groups, especially large ones.

## PRODUCT DESCRIPTION:

1. Massager
2. U-shaped head
3. Finger-shaped head
4. Ball-shaped head
5. Flat head
6. ON/OFF button
7. Intensity level LED (red)
8. Charge port
9. Charger
10. Charger LED

## CHARGING:

- Before using the device for the first time, charge it for three hours so the battery is fully charged. Connect the massager (1) to the power supply using the charger (9).

- The device takes approximately 2-3 hours to charge.

- When the massager (1) is charging, the charger LED (10) should be red.

- When the charger LED (10) turns green, the massager (1) is fully charged.

- The battery lasts for approximately 1-2 hours, depending on the speed chosen and the pressure applied.

## FUNCTIONS:

- The Gymform® Percussion Massager has six intensity levels.

## BENEFITS / TIME OF USE

Depending on the desired benefit, we recommend using the Gymform® Percussion Massager for a minimum of 15 seconds and a maximum of two minutes per muscle group. A full-body session should not last longer than 15 minutes. The device can be used 2 or 3 times on the same muscle group in the space of a single day.

### 1. Activating the Muscle

Before starting a training session or physical activity, use the device on the muscle for 30 seconds to instantly activate it.

### 2. Reactivating the Muscle

To reactivate the muscles during a training session or physical activity, use the device for 15 seconds on the muscle group. This will revitalize the muscles, alleviate muscle spasms and fatigue, and keep the muscles activated.

### 3. Pain Relief and Recovery

To reduce pain, use the Gymform® Percussion Massager for up to two minutes on the affected muscle group. We recommend using it for 30-second intervals with 5-10 seconds of rest between intervals, depending on the level of pain and sensitivity.

#### 4. Relaxation

To relax a particular muscle group, use the device for 30 seconds. This should be sufficient to relax the muscles in question.

#### Directions:

- Take the device out of the box.
- Attach the desired head, as shown in Fig. B.
- Press the ON/OFF button (6) to turn on the device (1). When the device is turned on, the intensity level LED (7) should indicate intensity level 1. Every time you press the button, the intensity level will increase and one of the six LEDs will turn on until the highest level (level 6) is reached. If you press the button again after reaching level 6, the device will switch off. Select the intensity level that best suits your requirements.
- To start the massage, place the head against the area you want to treat. Massage the area in accordance with the indications in the Benefits/Time of Use section above.
- Use the device only on soft tissue without causing pain or discomfort.
- Use the device only on dry, clean areas of skin. We recommend not massaging through clothing.
- Gently press and move the device over the skin. If you come across a painful area, keep the device focused on that spot.
- When you finish using the Gymform® Percussion Massager, hold down the ON/OFF button until the device switches off. When the device switches off, the intensity level LEDs (7) should also turn off.
- The device switches off automatically after 15 minutes.
- To remove the head, simply pull it off gently.

#### Cleaning Instructions:

Clean the device using a slightly damp cloth and then use a soft cloth to dry it. After cleaning, put the device away. Use the cardboard box the device came in to take it with you when you travel.

#### Specifications:

Input: 12 V 

Power: 30 W

#### Technical specifications of the adapter

Input: 100-240 V

Output: 12 V  1A

Frequency: 50-60 Hz

#### Disposal of materials

A symbol of a crossed-out wheeled bin means you should find out about and follow local regulations about disposing of this kind of product. Do not dispose of this product as you would other household waste. Dispose of this device in accordance with the corresponding local regulations. Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and/or human health and should be recycled properly.

#### Quality Warranty:

This product is covered by a warranty against manufacturing defects subject to the time periods stipulated by the legislation in force in each country.

This warranty does not cover damages resulting from inadequate use, negligent commercial use, abnormal wear and tear, accidents or improper handling.

#### Australia & New Zealand only

This product has the benefit of certain Consumer Guarantees. These are prescribed by the Australian Consumer Law & the New Zealand Consumer Guarantees Act 1993 both of which provide protection for consumers. There is no express warranty for this product in Australia or New Zealand. The above paragraph refers to other countries.

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ISL GmbH
Best Direct (International) Benelux B.V.
Best Direct International Ltd,
ISL ITALY S.r.l.