ENGLISH WARNINGS

BEFORE USING THE OR CHARGING THE Prime, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL, ON THE CHARGER, AND ON THE Prime. The Prime is intended for commercial or home use. We strive to make the Prime as safe as possible. This is an advanced mechanical device with electric components. If the Prime and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not place your finger or any objects near the metal piston above the attachment while the Prime is in use or immediately after use. The metal piston can become extremely hot. Do not use any attachments other than the ones provided by Theragun and only as instructed by Theragun. Do not use the Prime on your head or near your genitals. Do not use the Prime above the waist during pregnancy. Do not get hair caught in any moving part of the Prime. Do not use the Prime above your Adam's Apple or C4. Do not get the Prime wet. Only clean the Prime with a magic eraser sponge or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the Prime. Never store the Prime at temperatures above 40 C/104 F. Only charge the Prime when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the Prime or its charger to water. For long-term storage, the Prime should be fully charged. This product is not intended to diagnose, treat or prevent any diseases. When using the Prime, these basic precautions should always be followed including:

1. **USE ONLY AS INSTRUCTED.** Use the Prime as described in this Prime Instruction Manual. Use only Theragun recommended accessories and replacement parts. Do not carry out any maintenance other than what is shown in the Prime Instruction Manual or advised by Theragun.
2. **NOT FOR CHILDREN.** The Prime and charger are not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the Prime to be used as a toy. Close attention is necessary when used by or near children. Children should be supervised to ensure that they do not play with the Prime or charger.
3. **REDUCING RISK OF INJURY.** Charge the Prime only with the Theragun Prime charger. Other chargers may cause personal injury or damage. Do not connect the Prime to a different power supply plug or car cigarette lighter. The internal battery will be permanently disabled or damaged.
4. **CHARGING LOCATIONS.** The Prime should be charged indoors in a well ventilated, dry location. Do not charge the Prime outdoors, in a bathroom or within 10 feet (3.1 metres) of a bathtub or pool. Do not use the Prime or charger on wet surfaces and do not expose the charger to moisture, rain or snow. Do not use the Prime or charger in the presence of explosive atmospheres (gaseous fumes, dust or flammable materials). Sparks may be generated, possibly causing a fire.
5. **UNPLUGGING THE CHARGER.** Pull the plug, not the cable, to reduce the risk of damage to the electrical plug and cable. Never carry charger by its cable. Keep the cable from heat, oil and sharp edges. Do not use or charge with a damaged cable. Do

not stretch the charger cable or place the cable under strain. Keep cable away from heated surfaces. Do not handle the charger, including the charger terminals or Prime with wet hands.

6. **CHARGER CARE.** Unplug the charger when not in use. Unplug the charger from the socket when not in use for extended periods. Make sure the charger cable is located so that it will not be stepped on, tripped over, or otherwise subjected to damage or stress. Have a damaged charger replaced immediately.

7. **DO NOT OVERCHARGE.** Do not leave connect the Prime to the charger for more than an hour after the Prime has been fully charged. The Prime includes a protection system to avoid the risk of overcharging, however, overcharging the Prime may reduce its life over time. Only use with AC adapter: HXD302-1501500CI/ZZ, EA1023SAR.

8. **DO NOT BLOCK MOTOR VENTS.** Do not put any obstructions into vents. Do not use if any vent is blocked. Keep free of dust, lint, hair and anything that may reduce airflow. Keep hair, loose clothing, and fingers from vents and moving parts.

9. **DO NOT BURN OR INCINERATE THE Prime.** The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when battery is burned.

10. **DO NOT CRUSH, DROP OR DAMAGE THE Prime CHARGER.** Do not use a charger that has received a sharp blow, been dropped, run over, or damaged in any way.

11. **BATTERY CHEMICALS CAUSE SERIOUS BURNS.** Never allow the internal battery to come into contact with the skin, eyes or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water and rinse with vinegar. If eyes are exposed to battery chemicals, immediately flush with water for 20 minutes and seek medical attention. Remove and dispose of contaminated clothing.

12. **DO NOT OPERATE UNDER BLANKET AND PILLOW.** Excessive heating can occur and cause fire, electric shock, or injury to persons.

13. **STORING THE Prime AND CHARGER.** Store in a cool, dry place. Do not store the Prime where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.

14. **DISPOSING OF BATTERIES.** The internal Theragun Lithium-ion batteries for the Prime are more environmentally friendly than some other types of batteries (e.g. nickel-cadmium). Always dispose of the Prime according to government, regional and local regulations. Contact a recycling agency in your area for recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the Prime's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.

15. **DO NOT DISASSEMBLE.** Disassembly or incorrect reassembly may result in the risk of electric shock, fire or exposure to battery chemicals. The warranty will be void if the Prime or charger are disassembled or if any parts have been removed.

16. **SERVICE.** If the Prime or charger is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use and contact Theragun by emailing info@therabody.com for repair or service. Do not attempt to repair or disassemble the Theragun Prime which may result in an electric shock or fire.

17. **DO NOT USE WHILE BATHING OR IN THE SHOWER, BATHTUB OR SINK.** Do not place or store the Prime where it can fall or

be pulled into a bathtub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.

18. THERMAL LIMITER. The Prime employs an automatic reset thermal limiter that shuts off the Prime to prevent overheating and fire.

LIMITED WARRANTY

For full warranty information, please visit www.therabody.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address: Therabody - Warranty Attn: Customer Service 3003 W. Olympic Blvd., Suite 106-140 Los Angeles, CA 90006 Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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