

**ANT**  
**Flora Wave Digital Scale**  
 PERSONAL DIGITAL WEIGHING SCALE

**WARRANTY**

Ant Flora Wave Digital Scale comes with 30 Months warranty.

**Warranty Card**

Customer Name : \_\_\_\_\_  
 Product Name : \_\_\_\_\_ Model No.: \_\_\_\_\_  
 Address: \_\_\_\_\_ Date of Purchase: \_\_\_\_\_  
 Dealer Name : \_\_\_\_\_ Phone No.: \_\_\_\_\_



\*Notice: Any alteration of the clauses and contents and copies are invalid.  
 \*All product names, logos, and brands are property of their respective owners.  
 All company, product and service names used here are for identification purposes only.  
 Use of these names, logos, and brands does not imply endorsement.

**ANT**  
**Flora Wave Digital Scale**  
 PERSONAL DIGITAL WEIGHING SCALE

**USER MANUAL & WARRANTY CARD**



Your new electronic personal scale has been designed to accurately indicate your weight gain or loss over a period of time and should give many years of service under normal use.

**Specifications**

- Battery : 2x "AAA" Batteries (not included)
- Product Size : 26 x 26 cm
- Screen : 60 x 26 mm
- Pattern Type : Silk Screen
- Thickness of Glass : 5 mm
- Glass R-Angle : R2
- Material : Tempered Glass + ABS Plastic
- Sensor : High-Precision Strain Sensor
- Capacity : 180kg
- Graduation : 50g
- Mini Effective : 5kg
- Weight Mode : Automatic
- Auto Reset Zero : Yes
- Auto Switch Off : Yes
- Unit : kg
- Hint Of Out : o-Ld
- Hint Of Wrong : --
- Circuit Breaker Display : LCER
- Symbol Lo Boot Low-Voltage Detection, Low Voltage

**In The Box**

Weighing Scale x1, User Manual x1

**Key press**

- key functions: unit conversion switch  
 Program to load switch

**How to weight**

1. Remove the insulation film, put the scale on a hard & flat surface (avoid carpet or soft surface), wait 2 seconds before weighing.
2. Set up unit prior to use, unit can not be converted when weight.
3. Stand on the scale in balance, you can read the result automatically when the weight is set after 3 seconds and automatically switch off after 8 seconds.

Note: Before putting the scale in the flat & hard surface, the scale is automatically switch off, you must weight automatically first.

**ADVICE FOR USE AND CARE**

- When unused the scale should be left unloaded to avoid over assumption of power and failure of work.
- Try to stand in the middle place to avoid slip.
- Stand on the scale in good balance to avoid failure of lock and incorrect illustration.
- No heavy knock in the scale.
- Keep the scale dry. Do not use it in areas with water stains. Do not rinse with water.
- Handle with care, when using.

