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25.5cm			12cm
8.5cm	8.5cm	8.5cm	
<h2 style="margin: 0;">USER MANUAL</h2> <p style="margin: 0;">Smart USB Teeth Whitening Kit</p>	<p>Smart teeth whitening mouth tray led light adopts special cold light technology, combined with mouthpiece and high intensity cold light. Blue ray reacts with the professional teeth whitening gel, which can effectively remove the stains on the tooth surface.</p> <p>The soft mouthpiece is designed based on the principles of human oral engineering. You will feel more comfortable in the course of teeth whitening process. When you bite the mouth piece, it will completely fit and cover your teeth.</p> <p>Before starting your treatment, it is very important that you know how it works and what needs to be prepared. Please read through these instructions carefully before starting your treatment!</p>	<p>Start to Whiten</p> <div style="display: flex; flex-wrap: wrap; justify-content: space-around;"> <div style="text-align: center; margin: 5px;">  Rinse mouth with water </div> <div style="text-align: center; margin: 5px;">  Apply the gel on tooth surface </div> <div style="text-align: center; margin: 5px;">  Start by connecting to Apple/Android phone or USB port </div> <div style="text-align: center; margin: 5px;">  16 minute treatment </div> <div style="text-align: center; margin: 5px;">  Rinse mouth with water </div> <div style="text-align: center; margin: 5px;">  Rinse mouthpiece </div> </div>	

25.5cm			12cm
8.5cm	8.5cm	8.5cm	
<p>HOW TO USE</p> <ol style="list-style-type: none"> 1. Rinse mouth piece and your mouth before use. 2. Looking in the mirror or cell phone to compare your teeth with shade guide paper to record which stage you are in. 3. Twist teeth whitening pen, until the gel flows to the brush. Apply the gel evenly on your teeth, both upper and lower. Do not brush too much gel, more gel won't give you a faster or better result. Each time use about 0.5ml-0.7ml. Avoid applying the gel directly on your gum or inside of your lips! It will create a burning feeling. This is nothing to be alarmed about as the burning feeling will go away within 24 hours. Just make sure to rinse with water. 4. Make sure full charge before first use. Put mouth piece into your mouth and bite it, then turn on the light. Take deep breaths before treatment to decrease saliva production. (Excessive amount of saliva will dilute the gel's concentration.) 5. The whitening device will turn off automatically after 16 minutes. Take it out, rinse the mouth piece and your mouth with water and compare your teeth with shade guide paper again to see how many stages you have improved. 6. Do not eat or drink anything except water for 30 minutes after you finished the treatment. 	<p>INGREDIENTS</p> <p>Carbamide Peroxide, Glycerol, Propylene Glycol, Deionized Water, Carbomer, Carboxymethyl, Poly vinylpyrrolidone, Sodium Hydroxide, Menthol.</p> <p>RANGE OF APPLICATION</p> <p>Teeth Stains from coffee, tobacco, food coloring, aging and other factors that may cause discoloration.</p> <p>STORAGE</p> <p>Avoid sunlight and keep in cool and dry area. The gel pen can and should be kept in refrigerator.</p> <p>NOTES</p> <p>Avoid to apply the gel on the gum directly as possible as you can. It will create a burning feeling. This is nothing to be alarmed about as the burning feeling will go away within 24 hours. Just make sure to rinse with water.</p> <p>TERMS AND CONDITIONS</p> <p>Do not use if you are under the age of 18. Do not use if you have overly sensitive teeth. The product should only be used to whiten your teeth. Do not use if your teeth are decaying or</p>	<p>loose. Whitening product will not whiten any existing dental restoration such as fillings, crowns, veneers, bridges. Do not use if you are pregnant or lactating. Do not use if allergic to the ingredients listed above. Keep out of reach of children. Not to be swallowed.</p> <ol style="list-style-type: none"> 1. Avoid directly applying teeth whitening gel to gum area - wipe off excess. 2. If gel is on your face or clothes, you should wash them without any delay. 3. Please put cap on the top of pen after finishing teeth whitening treatment, which can avoid active ingredient volatilization. 4. Please store our product in a place which is out of children's reach. <p>CAUTION</p> <ol style="list-style-type: none"> 1. Not suitable for ceramic or false teeth. 2. Not suitable for teeth discoloration caused by lesion or medicine. 3. Not suitable for severe tetracycline and decayed teeth. 4. Not suitable for defective enamel, exposed dentin and damaged teeth. 5. Not suitable to children under 18 years old and pregnant women. 	