## ENGLISH WARNINGS

BEFORE USING OR CHARGING THE THERAGUN ELITE, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL, ON THE CHARGER, AND ON THE THERAGUN ELITE. The Theragun Elite is intended for commercial or home use. We strive to make the Theragun Elite as safe as possible. This is an advanced mechanical device with electric components. If the Theragun Elite and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not place your finger or any objects near the metal piston above the attachment while the Theragun Elite is in use or immediately after use. The metal piston can become extremely hot. Do not use any attachments other than the ones provided by Therabody and only as instructed by Therabody. Do not use the Theragun Elite on your head or near your genitals. Do not use the Theragun Elite above the waist during pregnancy. Do not get hair caught in any moving part of the Theragun Elite. Do not use the Theragun Elite above your Adam's apple or C4. Do not get the Theragun Elite wet. Only clean the Theragun Elite with a magic eraser sponge or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the Theragun Elite. Never store the Theragun Elite at temperatures above 40 C/104 F. Only charge the Theragun Elite when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the Theragun Elite or its charger to water. For long-term storage, the Theragun Elite should be fully charged. This product is not intended to diagnose, treat or prevent any diseases. When using the Theragun Elite, these basic precautions should always be followed including:

1. **USE ONLY AS INSTRUCTED.** Use the Theragun Elite as described in this Theragun Elite Instruction Manual. Use only Therabody recommended accessories and replacement parts. Do not carry out any maintenance other than what is shown in the Theragun Elite Instruction Manual or advised by Therabody.
2. **NOT FOR CHILDREN.** The Theragun Elite and charger are not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the Theragun Elite to be used as a toy. Close attention is necessary when used by or near children. Children should be supervised to ensure that they do not play with the Theragun Elite or charger.
3. **REDUCING RISK OF INJURY.** Charge the Theragun Elite only with the Theragun Elite charger. Other chargers may cause personal injury or damage. Do not connect the Theragun Elite to a different power supply plug or car cigarette lighter. The internal battery will be permanently disabled or damaged. Only use the appropriate Theragun Wireless Charger to wirelessly charge the Theragun Elite.
4. **CHARGING LOCATIONS.** The Theragun Elite should be charged indoors in a well ventilated, dry location. Do not charge the Theragun Elite outdoors, in a bathroom or within 10 feet (3.1 metres) of a bathtub or pool. Do not use the Theragun Elite or charger on wet surfaces and do not expose the charger to moisture, rain or snow. Do not use the Theragun Elite or charger in the presence of explosive atmospheres (gaseous fumes, dust or flammable materials). Sparks may be generated, possibly causing a fire.
5. **UNPLUGGING THE CHARGER.** Pull the plug, not the cable, to reduce the risk of damage to the electrical plug and cable. Never carry the charger by its cable. Keep the cable from heat, oil and sharp edges. Do not use or charge with a damaged cable. Do not stretch the charger cable or place the cable under strain. Keep cable away from heated surfaces. Do not handle the charger, including the charger terminals or Theragun Elite with wet hands.
6. **CHARGER CARE.** Unplug the charger when not in use. Unplug the charger from the socket when not in use for extended periods. Make sure the charger cable is located so that it will not be stepped on, tripped over, or otherwise subjected to damage or stress. Have a damaged charger replaced immediately.
7. **DO NOT OVERCHARGE.** Do not leave the Theragun Elite connected to the charger for more than an hour after the Theragun Elite has been fully charged. The Theragun Elite includes a protection system to avoid the risk of overcharging. However, overcharging the Theragun Elite may reduce its life over time. Only charge with AC Adapter : EA1046SAR or with the Theragun Elite Wireless Charger Stand.
8. **DO NOT BLOCK MOTOR VENTS.** Do not put any obstructions into vents. Do not use if any vent is blocked. Keep free of dust, lint, hair and anything that may reduce airflow. Keep hair, loose clothing, and fingers from vents and moving parts.

9. DO NOT BURN OR INCINERATE THE Theragun Elite. The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when battery is burned.
10. DO NOT CRUSH, DROP OR DAMAGE THE Theragun Elite CHARGER. Do not use a charger that has received a sharp blow, been dropped, run over, or damaged in any way.
11. BATTERY CHEMICALS CAUSE SERIOUS BURNS. Never allow the internal battery to come into contact with the skin, eyes or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water, then rinse with vinegar. If eyes are exposed to battery chemicals, immediately rinse with water for 20 minutes and seek medical attention. Remove and dispose of contaminated clothing.
12. DO NOT OPERATE UNDER BLANKET AND PILLOW. Excessive heating can occur and cause fire, electric shock, or injury to persons.
13. STORING THE Theragun Elite AND CHARGER. Store in a cool, dry place. Do not store the Theragun Elite where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.
14. DISPOSING OF BATTERIES. The internal Theragun Lithium-ion batteries for the Theragun Elite are more environmentally friendly than some other types of batteries (e.g. nickel-cadmium). Always dispose of the Theragun Elite according to government, regional and local regulations. Contact a recycling agency in your area for recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the Theragun Elite's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.
15. DO NOT DISASSEMBLE. Disassembly or incorrect reassembly may result in the risk of electric shock, fire or exposure to battery chemicals. The warranty will be void if the Theragun Elite or charger are disassembled or if any parts have been removed.
16. SERVICE. If the Theragun Elite or charger is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use and contact Therabody by emailing [info@therabody.com](mailto:info@therabody.com) for repair or service. Do not attempt to repair or disassemble the Theragun Elite which may result in an electric shock or fire.
17. DO NOT USE WHILE BATHING OR IN THE SHOWER, BATHTUB OR SINK. Do not place or store the Theragun Elite where it can fall or be pulled into a bathtub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.
18. THERMAL LIMITER. The Theragun Elite employs an automatic reset thermal limiter that shuts off the Theragun Elite to prevent overheating and fire.

This device can be used by children aged from 8 years and above, and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the device in a safe way and understand the hazards involved. Children shall not play with the device. Cleaning and user maintenance shall not be made by children without supervision.

The Theragun Elite is only to be used with the power supply unit provided with the device. Use only with Theragun Elite Adapter EA1046SAR, 45W. 

This device contains batteries that are only replaceable by Therabody authorized personnel.

This device contains batteries that are non-replaceable.

## LIMITED WARRANTY

For full warranty information, please visit [www.therabody.com/warranty](http://www.therabody.com/warranty). To request a copy of the warranty by mail, you may send a request to the following address: Therabody - Warranty Attn: Customer Service 3003 W. Olympic Blvd., Suite 106-140 Los Angeles, CA 90006. Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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