



PAUL HEWITT

PHREP tape measure

Suitable for all PAUL HEWITT PHREPS

Makes it easy to find the right size strap for you. The PHREP tape measure gives you the opportunity to find out the perfect size for your PAUL HEWITT PHREPS in four easy steps.




1. Print

Print DIN A4 in original size.
(Actual size or 100%)



2. Cut out

Cut out the tape measure and cut a slit at the point that has been marked.




3. Measure the length

Place the tape measure around your wrist and stick the point in the slit designed for that purpose. Make sure the piece of paper is not too tight against your wrist.



4. Read the number and place your order

Read the number that is located next to the opening. Now compare the number with the table below and order your desired PHREP at www.paul-hewitt.com.



Size table (valid for all PAUL HEWITT PHREPS)

Total length PHREPS	Circumference of your wrist
XS (16cm)	14-15 cm
S (17cm)	15-16 cm
M (18cm)	16-17 cm
L (19cm)	17-18 cm
XL (20cm)	18-19 cm
XXL (21cm)	19-20 cm
XXXL (22cm)	20-21 cm



Reference size

In order to check the precision of the tape measure, place a 1 Euro coin on top.



Make incision here